



Ultimate Perfection



DENTAL INSTRUMENTS CATALOGUE



TABLE OF CONTENTS

Page 01-05

Scissors

01

Page 06-14

Tweezers, Forceps & Needle Holder

02

Page 15-33

Extracting Forceps

03

Page 34-73

Mirror Handles, Scalars & Knives

04

Page 74-83

Root Elevators

05

Page 84-110

Implantology Instruments & Mouth Gags

06

07

Dental Syringes

Page 111 - 115

08

Crown, Ligature & Laboratory Instruments

Page 116 - 129

09

Articulators

Page 130 - 131

10

Orthodontic Instruments

Page 132 - 149

11

Impression Trays

Page 150 - 154

12

Hollow Ware

Page 155 - 164

13

Surgical Burs

Page 165 - 170




Page 01–05

- SCISSORS
- GUM & OPERATING SCISSORS
- SCISSORS WITH TUNGSTEN CARBIDE INSERTS
- MICRO SCISSORS



Gum & Operating Scissors



Scissors, One Blade Serrated

- One blade is with fine serrations to avoid tissue or ligature slippage.




Iris
11.5 cm



Goldman-Fox
13 cm



MM-SC-046
Goldman-Fox
13 cm



Kelly
16 cm



MM-SC-049
Locklin
16 cm

MM-SC-050
Locklin
16 cm

Super-Cut

- Cutting Blade horned to "razor" sharpness with opposing blade finished with fine serrations to avoid tissue slippage.



Iris
11.5 cm



MM-SC-053
La Grange
11.5 cm




Goldman-Fox
13 cm




Scissors With Tungsten Carbide Inserts




Iris
11.5 cm




Metzenbaum
14.5 cm




Kelly
16 cm




MM-SC-062
La Grange
11.5 cm




Goldman-Fox
13 cm




Strabismus
11.5 cm




MM-SC-067
15 cm




MM-SC-068
15 cm



MM-SC-069
18 cm



MM-SC-070
18 cm



MM-SC-071
18 cm



Micro Scissor
Stainless Steel



Micro Scissor
Stainless Steel




Micro Scissor
Stainless Steel
for ligature, one blade serrated


Micro Scissors with Titanium Stainless Steel Tips




MM-SC-072
18 cm




MM-SC-073
18 cm



MM-SC-074
18 cm



MM-SC-075
15.5 cm





MM-SC-076
15.5 cm







Page 06–14

- TWEEZERS
- DRESSING & DISSECTING FORCEPS
- HAEMOSTATIC FORCEPS
- NEEDLE HOLDERS
- NEEDLE HOLDERS WITH TUNGSTEN CARBIDE INSERTS
- MICRO NEEDLE HOLDERS & FORCEPS



Dissecting & Tissue Forceps








Tissue Forceps



Splinter Forceps

Towel Forceps



Haemostatic Forceps





Needle Holders



MM-NH-001
Derf
12.5 cm

MM-NH-002
Olsen-Hegar
14 cm

MMNH-003
Baumgartner
14.5 cm

MM-NH-004
Crilewood
15 cm


Mayo-Hegar
MM-NH-005 14 cm
MM-NH-006 15 cm
MM-NH-007 16 cm
MM-NH-008 18 cm
MM-NH-009 20 cm



Mathieu
MM-NH-010 14 cm
MM-NH-011 17 cm
MM-NH-012 20 cm



Mathieu-Kocher
MM-NH-013 14 cm
MM-NH-014 17 cm



MM-NH-015 14 cm
MM-NH-016 17 cm









MM-NH-017 14 cm
MM-NH-018 17 cm




Boynton
MM-NH-019
12 cm




Gillies with Scissor
MM-NH-020 Left handed
MM-NH-021 Right handed
16 cm




Needle Holders




Needle Holders with Tungsten Carbide inserts




Micro Needle Holders Stainless Steel



Micro Needle Holders Stainless Steel with Tungsten Carbide Inserts



Micro Forceps Stainless Steel





Page 15–33

- EXTRACTING FORCEPS




Fig. 1
MM-EF-001
Upper Centrals &
Canines / Cuspids



Fig. 2
MM-EF-002
Upper Laterals &
Canines / Cuspids



Fig. 3
MM-EF-003
Upper Incisors
& Canines



Fig. 4
MM-EF-004
Lower Incisors
& Cuspids



Fig. 7
MM-EF-005
Upper Premolars / Bicuspid




Fig. 8
MM-EF-006
Lower Premolars / Bicuspid



Fig. 13
MM-EF-007
Upper Molars, Right



Fig. 17
MM-EF-008
Upper Molars, Left




Fig. 18
MM-EF-009
Upper Molars



Fig. 18A
MM-EF-010
Upper Molars



Fig. 19
MM-EF-011
Upper Third Molars



Fig. 21
MM-EF-012
Lower Molars



Extraction Forceps Standard English Pattern



Fig. 22
MM-EF-013
Lower Molars




Fig. 23
MM-EF-014
Lower Molars
Right



Fig. 24
MM-EF-015
Lower Molars
Left




Fig. 29
MM-EF-016
Upper Roots




Fig. 29N
MM-EF-017
Upper Roots



Fig. 30
MM-EF-018
Upper Roots




Fig. 31
MM-EF-019
Upper Roots



Fig. 32
MM-EF-020
Lower Molars



Fig. 33
MM-EF-021
Lower Roots



Fig. 33A
MM-EF-022
Lower Roots



Fig. 33M
MM-EF-023
Lower Roots



Fig. 41
MM-EF-024
Upper Roots




Fig. 44
MM-EF-025
Upper Roots



Fig. 45
MM-EF-026
Lower Roots




Fig. 46L
MM-EF-027
Lower Roots




Fig. 51
MM-EF-028
Upper Roots




Fig. 51A
MM-EF-029
Upper Roots




Fig. 51M
MM-EF-030
Upper Roots



Fig. 52
MM-EF-031
Upper Roots



Fig. 52A
MM-EF-032
Upper Teeth
& Resection



Fig. 55
MM-EF-033
For Separating
Upper Molars



Fig. 56
MM-EF-034
For Separating
Lower Molars



Fig. 67
MM-EF-035
Upper Third
Molars




Fig. 67A
MM-EF-036
Upper Third
Molars

Extraction Forceps Standard English Pattern



Fig. 91
MM-EF-037
For Separating
Lower Molars



Fig. 99
MM-EF-038
Lower Molars



Fig. 99 1/4
MM-EF-039
Lower Molars



Fig. 99 1/2
MM-EF-040
Lower Molars




Fig. 145
MM-EF-041
Very Fine Roots
Lower



Fig. 151 1/2
MM-EF-042
Very Fine Roots
Upper





Fig. 349
MM-EF-043
Very Fine Roots
Upper




Fig. 359
MM-EF-044
Very Fine Roots
Lower



MD 1
MM-EF-045
Incisors, Bicuspid & Roots
Upper



MD 3
MM-EF-046
Incisors, Bicuspid & Roots
Lower



MD 2
MM-EF-047
Upper Molars



MD 4
MM-EF-048
Lower Molars



Fig. 68
MM-EF-049
Lower Roots




Fig. 68A
MM-EF-050
Lower Roots



Fig. 73
MM-EF-051
Lower Molars



Fig. 74
MM-EF-052
Lower Roots



Fig. 74N
MM-EF-053
Lower Roots



Fig. 75
MM-EF-054
Lower Premolars / Bicuspid



Fig. 79
MM-EF-055
Lower Third Molars



Fig. 79A
MM-EF-056
Lower Third Molars



Fig. 86A
MM-EF-057
Lower Molars



Fig. 86C
MM-EF-058
Lower Molars



Fig. 87
MM-EF-059
Lower Molars



Fig. 89
MM-EF-060
Right Upper Molars



Fig. 90
MM-EF-061
Left Upper Molars

Extraction Forceps With English Pattern




Fig. 1
MM-EF-001
Upper Centrals &
Canines / Cuspids




Fig. 7
MM-EF-005
Upper Premolars / Bicuspid



Fig. 13
MM-EF-007
Upper Molars, Right



Fig. 17
MM-EF-008
Upper Molars, Left




Fig. 18
MM-EF-009
Upper Molars



Fig. 22
MM-EF-013
Lower Molars



Fig. 33
MM-EF-021
Lower Roots



Fig. 51A
MM-EF-029
Upper Roots



Fig. 67A
MM-EF-036
Upper Third
Molars




Fig. 79
MM-EF-043
Lower Third Molars



MM-EF-062
Set of 10 pieces



Fig.3
MM-EF-063
Lower Premolars




Fig.5
MM-EF-064
Lower Roots




Fig.6
MM-EF-065
Lower Roots



Fig.7
MM-EF-066
Upper Incisors




Fig.51S
MM-EF-067
Lower Incisors and Canines




Fig.137
MM-EF-068
Upper Molars, left




Fig.139
MM-EF-069
Upper Molars, Right



MM-EF-070
Set of 7 pieces

Extraction Forceps With Anatomical Grip




Fig. S-31
MM-EF-071
Upper Premolars




Fig.34
MM-EF-072
Upper Incisors & Cuspids




Fig.35
MM-EF-073
Upper Incisors & Cuspids,
Deep Gripping




Fig. 79N
MM-EF-074
Lower Wisdom Teeth




Fig. 46N
MM-EF-075
Lower Premolars



Fig.36
MM-EF-076
Lower roots
Incisors & Cuspids



Fig.36N
MM-EF-077
Lower roots
Incisors & Cuspids



Fig.59
MM-EF-078
Lower Root Fragments






Fig. 80
MM-EF-079
Felsch
Upper Wisdom Teeth
with Retaining Device




Fig. 67L
MM-EF-080
Upper Molars
Left



Fig. 67R
MM-EF-081
Upper Molars
Right




Fig. 67N
MM-EF-082
Upper Molars




Fig. 51B
MM-EF-083
Upper Molars




Fig. 97
MM-EF-084
Upper Root Fragments




Fig. 97B
MM-EF-085
Upper Root Fragments



Fig. 122 L
MM-EF-086
Lower Molars & Wisdom Teeth
Left



Fig. 122 R
MM-EF-087
Lower Molars & Wisdom Teeth
Right




Fig. 133 L
MM-EF-088
Lower Very Fine Roots
Left




Fig. 133 R
MM-EF-089
Lower Very Fine Roots
Right



Extraction Forceps "American Pattern"




Standard #1
MM-EF-090
For upper incisors & cuspids



Henahan #1 A
MM-EF-091
For upper incisors




#2
MM-EF-092
Upper & Lower Roots




Woodward #3 FS
MM-EF-093
For lower molars



#4
MM-EF-094
Upper & Lower Premolars




Physick #5
MM-EF-095
Lower third molars




Spaltzange #6
MM-EF-096
Splitting & Separating Forceps




#7
MM-EF-097
Upper Roots



MM-EF-098 # 8
Upper & Lower third molars



MM-EF-099 # 10 H
For upper Molars, Hook Handle




MM-EF-100 # 10 S
For upper Molars




MM-EF-101 # 13
Upper Incisors




MM-EF-102 # 14
Lower Incisors



MM-EF-103 # 15
For 1st & 2nd Lower Molars




"Cow Horn"
MM-EF-104 # 16
for 1st & 2nd Lower Molars




"Cow Horn"
MM-EF-105 # 16 SK
for children, for Lower Molars


Extraction Forceps "American Pattern"



MM-EF-106 # 17
MM-EF-107 # 17 Gold Plated
for 1st & 2nd Lower Molars



MM-EF-108 # 18 L
for 1st & 2nd Upper Molars, left




MM-EF-109 # 18 R
for 1st & 2nd Upper Molars, right




MM-EF-110 # 20
Upper Molars left




MM-EF-111 # 20 R
Upper Molars right



"Cow Horn"
MM-EF-112 # 23
for Lower Molars



MM-EF-113 # 24
for Upper Molars




MM-EF-114 # 27
for Upper Molars




MM-EF-115 # 29 S
Lower teeth for children



MM-EF-116 # 30 S
Upper teeth for children




MM-EF-117 # 32
for Upper Bicuspids



Parmly
MM-EF-118 # 32 A
for upper Bicuspids & Roots



MM-EF-119 # 34
Lower Anteriors & Roots



MM-EF-120 # 35
Upper Roots



MM-EF-121 # 36
Lower Anteriors & Roots




MM-EF-122 # 37
Lower Anteriors & Roots




Extraction Forceps "American Pattern"




MM-EF-123 # 38
Upper Anteriors & Roots




MM-EF-124 # 53 L
MM-EF-125 # 53 L Gold Plated
for upper Molars, left




MM-EF-126 # 53 R
MM-EF-127 # 53 R Gold Plated
for upper Molars, right




MM-EF 128 # 62
for upper & Lower Incisors




MM-EF-129 # 65
MM-EF-130 # 65 Gold Plated
for upper Incisors, Roots & Fragments




MM-EF-131 # 69
MM-EF-132 # 69 Gold Plated
for upper & Lower Fragments & Roots




MM-EF-133 # 85 A
for Lower Molars & Bicuspids




Nevius
MM-EF-134 # 88 L
for Upper Molars, left




Nevius
MM-EF-135 # 88 R
For upper Molars, right




MM-EF-136 # 89
For upper Molars, right




MM-EF-137 # 90
For upper Molars, left




Kells
MM-EF-138 # 99 A
Upper Incisors, Canines, Premolars




Kells
MM-EF-139 # 99 C
For upper Incisors & Bicuspid



Hull
MM-EF-140 # 01
Universal for Bicuspid & Deciduous




MM-EF-141 # 103
for Lower Incisors, Bicuspid & Roots




MM-EF-142 # 104
Upper third molars


Extraction Forceps "American Pattern"




Nevius
MM-EF-143 # 110 S
for Upper Bicuspid & Cuspid




Cryer
MM-EF-144 # 150
MM-EF-145 # 150 Gold Plated
for Upper Incisors, Bicuspid & Roots




Cryer
MM-EF-146 # 150 A
For upper Incisors, Cuspid & Bicuspid




Cryer
MM-EF-147 # 150 AS
For upper Incisors, Cuspid
serrated jaw




Cryer
MM-EF-148 # 150 XAS
For children, for Upper Incisors &
Cuspid, serrated jaw




Cryer
MM-EF-149 # 150 S
MM-EF-150 # 150 S Gold Plated
For children, for Upper Incisors, Bicuspid & Roots




Cryer
MM-EF-151 # 151
MM-EF-152 # 151 Gold Plated
For Lower Incisors, Bicuspid & Roots




Cryer
MM-EF-153 # 151 A
For Lower Incisors, Cuspid
& Bicuspid



Cryer
MM-EF-154 # 151 AS
For Lower Incisors & Cuspids, serrated jaw



MM-EF-155 # 151 XAS
For children, for Lower Incisors & Cuspids, serrated jaw




Cryer
MM-EF-156 # 151 S
MM-EF-157 # 151 S Gold Plated
For children, for Lower incisors, Bicuspids & Roots




MM-EF-158 # 201
Upper Incisors, Canines, premolars



MM-EF-159 # 203
Lower Incisors & Roots



MM-EF-160 # 210 H
For Upper Wisdoms, Hook Handle




MM-EF-161 # 210 S
For Upper Wisdoms




MM-EF-162 # 213
Upper Centrals & Canines




Extraction Forceps "American Pattern"



MM-EF-163 # 217
For Lower Molars




MM-EF-164 # 222
For Lower Wisdoms




MM-EF-165 # 286
For Upper Incisors, Bicusps & Roots




MM-EF-166 # 288
Lower Canines & Premolars




Cohen
MM-EF-167 # 310
Universal Forceps



Nevius
MM-EF-168 #320/10S
Upper Canines



Henahan
MM-EF-169 # 315/1A
Upper Incisors




Nevius
MM-EF-170 #325/10H
Upper Canines




Standard
MM-EF-171 # 330/4
Upper Third Molars




Woodward
MM-EF-172 # 335/FH
Lower Molars




MM-EF-173
17 SK
For children for Lower Wisdoms



MM-EF-174
23 SK
For children for Lower Wisdoms



MM-EF-175
150 SK
For Children, for Upper Incisors,
Bicuspids & Roots




MM-EF-176
151 SK
For Children, for Lower Incisors,
Bicuspids & Roots


Page

34–73



- MOUTH MIRROR & MIRROR HANDLES
- SCALPEL HANDLES & BLADES
- EXPLORERS SINGLE & DOUBLE ENDED
- PROBES SINGLE ENDED
- GRACEY CURETTES
- PERIODONTAL INSTRUMENTS
- SCALERS
- ENDODONTIC INSTRUMENTS
- FILLING INSTRUMENTS
- CAVITY PREPRATION TIPS & AMALGAM INSTRUMENTS
- AMALGAM PLUGGERS
- BURNISHERS
- COMPOSITE / PLASTIC FILLING INSTRUMENTS
- EXCAVATORS
- CARVERS / FINISHING KNIVES & CLEOID - DISCOID
- WAX MODELLING CARVERS
- CEMENT SPATULAS
- PLASTER KNIVES - ALGINATE & PLASTER SPATULAS
- WAX KNIVES



- Rhodium coated front surface mirror
- Wide stem offers comfortable cheek retraction
- Non glare satin finish



Pack of 12				
Fig.	Ø mm	Front surface		Rear surface
		Chrome	Rhodium	Magnifying
2	18	-	-	MM-MM-009
3	20	-	-	MM-MM-010
4	22	MM-MM-005	MM-MM-007	MM-MM-011
5	24	MM-MM-006	MM-MM-008	MM-MM-012
				Plane
				MM-MM-013
				MM-MM-014
				MM-MM-015
				MM-MM-016



Mirror handles solid



Mirror handles hollow





Fig. 1 MM-MH-001 Fig. 2 MM-MH-002 Fig. 3 MM-MH-003 Fig. 4 MM-MH-004

Fig. 1 MM-MH-005 Fig. 2 MM-MH-006 Fig. 3 MM-MH-007 Fig. 4 MM-MH-008 Fig. 5 MM-MH-009 Fig. 6 MM-MH-010 Fig. 7 MM-MH-011




Handles Selection For Scaler

Double Ended



Double Ended

Single Ended



Sterile
Individually Packed
Box of 100 Pieces




Fig. 10 MM-SB-001 Fig. 11 MM-SB-002 Fig. 12 MM-SB-003 Fig. 12D MM-SB-004 Fig. 13 MM-SB-005 Fig. 15 MM-SB-006 Fig. 15C MM-SB-007 Fig. 20 MM-SB-008 Fig. 22 MM-SB-009 Fig. 23 MM-SB-010 Fig. 24 MM-SB-011 Fig. 25 MM-SB-012





Fig. 3 MM-SH-001 Fig. 3G MM-SH-002 Fig. 4 MM-SH-003 Fig. 3L MM-SH-004 Fig. 7K MM-SH-005 Fig. 7 MM-SH-006 MM-SH-007 35 mm MM-SH-008 40 mm MM-SH-009 45 mm MM-SH-010 50 mm Fig. 5S MM-SH-011 Fig. 5A MM-SH-012 Fig. 5 MM-SH-013

Micro Scalpel Handles

Micro scalpel & micro mirror handles

Safety Model




2.5 cm

MM-SH-019

Scalpel Handle double Blade
15.5 cm

Scalpel Blade removing Forceps



1. Hold blade as shown.

2. Lift blade.

3. Push forward to remove blade.

Scalpel Handle Adjustable
Locks at Angles: 30°, 60°, 90°, 120°, 150°, 180°
MM-SH-023

Explorers double ended rigid



Fig. 23/17 **MM-EX-001**



Fig. 23/17A **MM-EX-002**

Explorers double ended rigid



Fig. 23/17 **MM-EX-003**



Fig. 23/17A **MM-EX-004**



Fig. 23/16 **MM-EX-005**

Explorers double ended Flexible



Fig. 2 **MM-EX-006**



Fig. 3 **MM-EX-007**



Fig. 5 **MM-EX-008**




Fig. TU17/23 **MM-EX-009**




Fig. 6/23 **MM-EX-010**

Explorers Double Ended

Expro (explorer-probe combination)





Probes double ended



Explorers single ended rigid

Explorers single ended flexible





Probes Single Ended

Explorers single ended rigid

Explorers single ended rigid



Fig. 2
MM-PR-024



Fig. 6
MM-PR-025



Fig. 8
MM-PR-026



Fig. 9
MM-PR-027



Fig. 10
MM-PR-028



Fig. 11
MM-PR-029



Fig. 12
MM-PR-030



Fig. 16
MM-PR-031



Fig. 17
MM-PR-032



Fig. 23
MM-PR-033



Fig. 1
MM-PR-034



Fig. 2
MM-PR-035



Fig. 3
MM-PR-036



Fig. 4
MM-PR-037



Fig. 5
MM-PR-038



Fig. 6
MM-PR-039



Fig. 7
MM-PR-040




Fig. 8
MM-PR-041



Fig. 9
MM-PR-042


Color Coded






Endodontic Instruments

Endodontic explorers





Solid handle



Endodontic excavators

Solid handle






Fig. 1
MM-EI-013

Perforation Repair Instruments




Fig. 1
MM-EI-014




Fig. 2
MM-EI-015




Fig. 1
MM-EI-016




Fig. 2
MM-EI-017

Hand Spreaders



Fig. 3
MM-EI-018

Solid handle




Fig. 1
MM-EI-019




Fig. 2
MM-EI-020




Fig. 4
MM-EI-021



Fig. 5
MM-EI-022




Fig. 6
MM-EI-023




Fig. 7
MM-EI-024




Fig. 1
MM-EI-025




Fig. 2
MM-EI-026



Fig. 3
MM-EI-027




Fig. 4
MM-EI-028

Endodontic Instruments



Fig.2/4
MM-EI-029
Molt




Fig.1
MM-EI-030
Ruddle
Left



Fig.2
MM-EI-031
Ruddle
Right




Fig.1
MM-EI-032
Ruddle
Left





Fig.2
MM-EI-033
Ruddle
Right

Endodontic Forceps

Micro forceps for removing broken files



Silver points and broken files removing forceps



Endodontic Tweezers






Fig. 1
MM-EI-041

Solid handle




Fig. 2
MM-EI-042




Fig. 3
MM-EI-043




Fig. 4
MM-EI-044




Fig. 1
MM-EI-045
Luks




Fig. 2
MM-EI-046
Luks




Fig. 3
MM-EI-047
Luks




Fig. 4
MM-EI-048
Luks




Fig. 1
MM-EI-049




Fig. 2
MM-EI-050




Fig. 3
MM-EI-051




Fig. 4
MM-EI-052

Excavators




Fig. 0 **MM-EC-001**




Fig. 1 **MM-EC-002**




Fig. 2 **MM-EC-003**




Fig. 3 **MM-EC-004**




Fig. 4 **MM-EC-005**



Solid handle



Fig. 6 **MM-EC-006**



Fig. 14 **MM-EC-007**



Fig. 17 **MM-EC-008**



Fig. 18 **MM-EC-009**



Fig. 19 **MM-EC-010**



Fig. 17W **MM-EC-011**



Fig. 18W **MM-EC-012**



Fig. 19W **MM-EC-013**




Fig. 0 **MM-EC-014**



Fig. 1 **MM-EC-015**



Fig. 2 **MM-EC-016**




Fig. 3 **MM-EC-017**



Fig. 4 **MM-EC-018**



Solid handle

Palma

Type A




Fig. 1 MM-EC-019



Fig. 2 MM-EC-020



Fig. 3 MM-EC-021



Fig. 4 MM-EC-022

Type C



Fig. 1 MM-EC-023




Fig. 2 MM-EC-024



Fig. 3 MM-EC-025



Fig. 4 MM-EC-026



Type B

Palma




Fig. 1 MM-EC-027



Fig. 2 MM-EC-028




Fig. 3 MM-EC-029




Fig. 4 MM-EC-030



Solid handle

Blacks




Fig. 77/78 MM-EC-031




Fig. 79/80 MM-EC-032

Ferrier




Fig. 26 MM-EC-033




Fig. 27 MM-EC-034



Fig. 28 MM-EC-035



Fig. 29 MM-EC-036



Burnishers



Acron / Wescott




Fig. 21 **MM-BR-001**




Fig. 21B **MM-BR-002**

White



Fig. 27/29 **MM-BR-003**

Freedman



Fig. 3 **MM-BR-004**

Solid Handle



Palma




Fig. 1 **MM-BR-011**




Fig. 2 **MM-BR-012**




Fig. 3 **MM-BR-013**




Fig. 1 **MM-BR-014**



Fig. 1L **MM-BR-015**




Fig. 18 **MM-BR-016**




Fig. 1 **MM-BR-017**




Fig. 2 **MM-BR-018**




Fig. 3 **MM-BR-019**




Fig. 4 **MM-BR-020**




Fig. 1 **MM-BR-021**




Fig. 2 **MM-BR-022**







Fig. 3 **MM-BR-023**









Fig. 4 **MM-BR-024**


Solid handle



Beavertail / White



White









Hollow handle

Acron / Wescott




White







Amalgam Pluggers

Solid handle






Markley



Solid handle





Black



Palma

Palma



Cavity Preparation Tips & Amalgam Instruments



MM-AC-001 Mini/Regular 20 cm




MM-AC-002 Regular/Large 20 cm




MM-AC-003 Regular/Jumbo 20 cm




MM-AC-004 Large/Jumbo 20 cm




MM-AC-005 Mini 15 cm



MM-AC-006 Regular 15 cm



MM-AC-007 Large 15 cm



MM-AC-008 Jumbo 15 cm



Fig. 10 **MM-AC-009** 16.5 cm
English pattern



Fig. 1 **MM-AC-010** 18.5 cm
German pattern



Fig. 1 **MM-AC-011** 15 cm
Swiss pattern



Fig. 1 **MM-AC-012** Regular 15 cm




Fig. 2 **MM-AC-013** Regular 15 cm



Fig. 3 **MM-AC-014** Regular 14 cm

For micro filling



MM-AC-015 Str
MM-AC-016 Cvd
13 cm

Amalgam well



MM-AC-017



Pluggers



Fig. 5A MM-CI-001



Fig. 3 MM-CI-002 Hollenback



Fig. 1 MM-CI-003 Super



Fig. OT MM-CI-004 Tanner



Fig. 0/1 MM-CI-005



Fig. ½ MM-CI-006

Burnishers



Fig. 21B MM-CI-007 Acom/Wescott



Fig. 1 MM-CI-008 Contact Former



Fig. 2 MM-CI-009 Contact Former



Fig. 2 MM-CI-010 Ladmore



Fig. 3 MM-CI-011 Ladmore



Fig. 27/29 MM-CI-012



Fig. 1 MM-CI-013



Fig. 1L MM-CI-014



Fig. 155 MM-CI-015

Composite Instruments



Fig. 156(A6) **MM-CI-016**



Fig. 8A **MM-CI-017**



Fig. 1 **MM-CI-018** Boghosian



Fig. 2 **MM-CI-019** Boghosian



Fig. 1 **MM-CI-020**



Fig. 2 **MM-CI-021**



Fig. 1 **MM-CI-022** Goldstein



Fig. 2 **MM-CI-023** Goldstein



Fig. 3 **MM-CI-024** Goldstein



Fig. 4 **MM-CI-025** Goldstein



Fig. 1 **MM-CI-026** Goldstein-Mini



Fig. 3 **MM-CI-027** Goldstein-Mini



Fig. 4 **MM-CI-028** Goldstein-Mini



Fig. 1 **MM-CI-029** Goldstein-Micro



Fig. 4/5 **MM-CI-030** Gregg



Fig. 0 **MM-CI-031** Heidemann



Fig. 1 **MM-CI-032** Heidemann



Fig. 2 **MM-CI-033** Heidemann



Fig. W3 **MM-CI-034**



Fig. 1 **MM-CI-035**



Fig. 6 **MM-CI-036**



Fig. 178 **MM-CI-037**



Fig. 49 **MM-CI-038**



Fig. 179 **MM-CI-039**




Fig. 6 MM-PF-001



Fig. 178 MM-PF-002



Fig. 179 MM-PF-003



Fig. 5/7 MM-PF-004

Gregg




Fig. 4/5 MM-PF-005

Woodson



Fig. 1 MM-PF-006



Fig. 2 MM-PF-007



Fig. 3 MM-PF-008

Heidemann




Fig. 0 MM-PF-009




Fig. 1 MM-PF-010




Fig. 2 MM-PF-011




Fig. 0 MM-PF-012



Fig. 1 MM-PF-013



Fig. 2 MM-PF-014



Fig. 3 MM-PF-015



Fig. 4 MM-PF-016




Fig. 5 MM-PF-017




Fig. 6 MM-PF-018




Fig. 7 MM-PF-019



Fig. 8 MM-PF-020




Fig. 9 MM-PF-021




Fig. 10 MM-PF-022













Fig. 11 MM-PF-023



Fig. 12 MM-PF-024





Plastic Filling Instruments


Hollow handle





Heidemann

Solid handle






Heidemann




Gregg


Applicators for calcium hydroxide



Applicators for medicaments



Scalers



Solid handle




Fig. 204
MM-SC-014




Fig. CI 2/3
MM-SC-021




Fig. 204S
MM-SC-015




Fig. 6
MM-SC-022
Crane-Kaplan




Fig. 204SD
MM-SC-016




Fig. 6/7
MM-SC-023
Hygienist




Fig. 31/32
MM-SC-017
Jacquette




Fig. UI5/33
MM-SC-024
Jacquette/Towner




Fig. 30/33
MM-SC-018
Jacquette




Fig. 33/H5
MM-SC-019
Jacquette




Fig. 2/3
MM-SC-025
Taylor




Fig. UI5/30
MM-SC-020
Jacquette/Towner




Fig. 7/8
MM-SC-026
Younger-Good

Scalers




Fig. 204
MM-SC-027




Fig. CI 2/3
MM-SC-034




Fig. 204S
MM-SC-028




Fig. 6
MM-SC-035
Crane-Kaplan




Fig. 204SD
MM-SC-029




Fig. 6/7
MM-SC-036
Hygienist




Fig. 31/32
MM-SC-030
Jacquette




Fig. U15/33
MM-SC-037
Jacquette/Towner




Fig. 30/33
MM-SC-031
Jacquette




Fig. 33/H5
MM-SC-032
Jacquette




Fig. 2/3
MM-SC-038
Taylor




Fig. U15/30
MM-SC-033
Jacquette/Towner





Fig. 7/8
MM-SC-039
Younger-Good


Solid Handle




Darby Perry




Standard











Standard





Palma







Hygienist



Jacquette




Towner




Hollow Handle


Hygienist









Jacquette




Standard



Standard



Towner



Periodontal chisels



Fig. 1 MM-PI-001 TG




Fig. 2 MM-PI-002 TG

Universal Curettes Columbia

Solid handle




Hollow handle



Solid handle

Hollow handle




Universal

Fig. 13R/14L

MM-PI-003

MM-PI-006




Anterior teeth

Fig. 2R/2L

MM-PI-004

MM-PI-007



Posterior teeth

Fig. 4R/4L

MM-PI-005

MM-PI-008

Hoe scalars


Solid handle



Hollow handle



Goldman-Fox



Solid handle

Hollow handle

Fig. 5

MM-PI-009

MM-PI-011




Fig. 6

MM-PI-010

MM-PI-012

Solid handle



Hollow handle



Anterior teeth




Fig. 1/2

Anterior teeth




Fig. 3/4

Anterior and
Premolars




Fig. 5/6

Premolars and
Molars




Fig. 7/8

Molars buccal and
lingual surfaces




Fig. 9/10

Posterior mesial
surfaces



Fig. 11/12

Posterior distal
surfaces




Fig. 13/14

Solid handle

	Solid handle
Fig. 1/2	MM-GC-001
Fig. 3/4	MM-GC-002
Fig. 5/6	MM-GC-003
Fig. 7/8	MM-GC-004
Fig. 9/10	MM-GC-005
Fig. 11/12	MM-GC-006
Fig. 13/14	MM-GC-007

MM-GC-008

Gracey set of 7 pieces

Hollow handle

	Hollow handle
Fig. 1/2	MM-GC-009
Fig. 3/4	MM-SC-010
Fig. 5/6	MM-GC-011
Fig. 7/8	MM-GC-012
Fig. 9/10	MM-GC-013
Fig. 11/12	MM-GC-014
Fig. 13/14	MM-GC-015

MM-GC-016

Gracey set of 7 pieces

Gracey Curettes



Fig. 1/2
MM-GC-017
Grey
Anterior teeth



Fig. 3/4
MM-GC-018
Black
Anterior teeth



Fig. 5/6
MM-GC-019
Yellow
Anterior and Premolars



Fig. 7/8
MM-GC-020
Green
Premolars and Molars



Fig. 9/10
MM-GC-021
Red
Molars buccal and lingual surfaces



Fig. 11/12
MM-GC-022
Orange
Posterior mesial surfaces



Fig. 13/14
MM-GC-023
Blue
Posterior distal surfaces

MM-GC-024
Set of 7 pieces
Fig. 1/2, 3/4, 5/6, 7/8, 9/10, 11/12, 13/14



Fig. 15/16
MM-GC-025
Coffee
Posterior mesial surfaces




Fig. 17/18
MM-GC-026
Light Blue
Posterior distal surfaces




Fig. 11/14
MM-GC-027
Blue
Mesial/Distal surfaces



Fig. 12/13
MM-GC-028
Orange
Mesial/Distal surfaces

- Ergonomic, large diameter silicone handle for better grip
- Ultra light weight
- Special grade stainless steel tips ensure longer life
- Rigid pattern, thicker tips than standard Gracey



MM-GC-043
Set of 7 pieces / Rigid

MM-GC-044
Set of 7 pieces / Standard




Fig. 3 **MM-CR-001** Hollenback




Fig. 2/3 **MM-CR-005** Frahm




Fig. 1 **MM-CR-002** Interproximal




Fig. 1 **MM-CR-006** Wards




Fig. 1 **MM-CR-003** Interproximal




Fig. 2 **MM-CR-007** Wards




Fig. 2 **MM-CR-004** Interproximal




Fig. 8 **MM-CR-008** Wiland

Hollow handle



Fig. 6T **MM-FK-001** Tanner



Fig. 2 **MM-FK-002** USC



Fig. 2S **MM-FK-003** USC

Cleoid - Discoid



Solid handle



Fig. 1 **MM-CD-001**



Fig. 3 **MM-CD-003**



Fig. 2 **MM-CD-002**



Fig. 4 **MM-CD-004**

Cleoid - Discoid




Fig. 4/5 **MM-CD-005**




Fig. 89/92 **MM-CD-007**



Fig. 3/6 **MM-CD-006**




Fig. 90B **MM-CD-008**

Wax Modelling Carvers

Hollow Handle



Fig. 1 MM-MC-001




Fig. 2 MM-MC-002



Fig. 4 MM-MC-004



Fig. 3 MM-MC-003




Fig. 5 MM-MC-005

Solid Handle



Fig. 1 MM-MC-006



Fig. 2 MM-MC-007



MM-MC-011
P.k. Thomas
Set of 5 Pieces



Fig. 4 MM-MC-009



Fig. 3 MM-MC-008



Fig. 5 MM-MC-010

Hollow Handle



Fig. 7 MM-MC-012 Beale



Fig. 31 MM-MC-013 Gritman



Fig. 10 MM-MC-015 Lecron



Fig. 10 MM-MC-014 Hylin



Fig. 10 MM-MC-016 Zahle

Solid Handle



Fig. 1 MM-MC-017 Vehe



Fig. 31 MM-MC-018 Gritman



Fig. 10 MM-MC-020



Fig. 12 MM-MC-019



Fig. 1 MM-MC-021 Beale



Fig. 1 MM-MC-022 Roach



MM-MC-023 MM-MC-024 MM-MC-025 MM-MC-026 MM-MC-027 MM-MC-028 MM-MC-029 MM-MC-030 MM-MC-031 MM-MC-032 MM-MC-033 MM-MC-034 MM-MC-035



MM-MC-036 MM-MC-037 MM-MC-038 MM-MC-039 MM-MC-040 MM-MC-041 MM-MC-042 MM-MC-043 MM-MC-044 MM-MC-045 MM-MC-046 MM-MC-047 MM-MC-048



Fig. 1 MM-CP-001



Fig. 2 MM-CP-002



Fig. 3 MM-CP-003



Fig. 1 MM-CP-004



Fig. 2 MM-CP-005



Fig. 3 MM-CP-006




Fig. 00 MM-CP-007




Fig. 0 MM-CP-008




Fig. 1 MM-CP-009




Fig. 2 MM-CP-010




Fig. 3 MM-CP-011



Fig. 4 MM-CP-012



Weston MM-CP-013



MM-CP-014

Plaster knives with wooden handle



MM-WK-012 Fahnestock 13 cm



MM-WK-014 Lessmann 13 cm



MM-WK-013 Fahnestock 17.5 cm



MM-WK-015 Lessmann 17.5 cm

Wax knives with wooden handle



MM-WK-016 Fahnestock 13 cm



MM-WK-018 Lessmann 13 cm



MM-WK-017 Fahnestock 17.5 cm



MM-WK-019 Lessmann 17.5 cm



MM-WK-020 Con scodellino 13 cm



MM-WK-021 Con scodellino 17.5 cm


Wax knives Ermert



MM-WK-022 Wooden handle



MM-WK-023 Plastic handle



MM-WK-024 Metal handle

Wax knives with metal handle



MM-WK-025 Fahnestock 13 cm



MM-WK-026 Fahnestock 17.5 cm



MM-WK-027 Lessmann 13 cm




MM-WK-028 Lessmann 17.5 cm

Plaster Spatulas

Flexible Spatulas



Rigid Spatulas



Page
74–83

○ ROOT ELEVATORS



Root Elevators

Luxators



Fig. 1
MM-RE-001



Fig. 2
MM-RE-002



Fig. 3
MM-RE-003



Fig. 4
MM-RE-004




Fig. 5
MM-RE-005



Fig. 6
MM-RE-006




Fig. 7
MM-RE-007



MM-RE-008
Set of 7

Bein




Fig. 1
MM-RE-009




Fig. 2
MM-RE-010



Fig. 3
MM-RE-011



Fig. 4
MM-RE-012



Fig. 5
MM-RE-013



Fig. 6
MM-RE-014

Apical




Fig. 301
MM-RE-015



Fig. 302
MM-RE-016
Left





Fig. 303
MM-RE-017
Right




Fig. 304
MM-RE-018

Seldin




Cryer



Coupland



Schmeckebier





Root Elevators Bien



Fig. 1
MM-RE-036



Fig. 2
MM-RE-037



Fig. 3
MM-RE-038



Fig. 4
MM-RE-039



Fig. 5
MM-RE-040



Fig. 6
MM-RE-041



Fig. 7
MM-RE-042



Fig. 8
MM-RE-043



Fig. 9
MM-RE-044



Fig. 10
MM-RE-045



Fig. 1
MM-RE-046



Fig. 2
MMRE-047



Fig. 3
MM-RE-048



Fig. 1
MM-RE-049



Fig. 2
MM-RE-050



Fig. 3
MM-RE-051



Fig. 1
MM-RE-052
Bayonet Distal



Fig. 2
MM-RE-053
Bayonet Distal




Fig. 3
MM-RE-054
Bayonet Distal



Fig. 1
MM-RE-055
Bayonet Mesial




Fig. 2
MM-RE-056
Bayonet Mesial




Fig. 3
MM-RE-057
Bayonet Mesial



Flohr



Fig. 1
MM-RE-058




Fig. 2
MM-RE-059
Left




Fig. 3
MM-RE-060
Right

Friedman



Fig. 1
MM-RE-061




Fig. 2
MM-RE-062
Left



Fig. 3
MM-RE-063
Right

Lindo-Levian



Fig. LLL
MM-RE-064




Fig. LLM
MM-RE-065




Fig. LLS
MM-RE-066




Fig. LLL
MM-RE-067




Fig. LLM
MM-RE-068




Fig. LLS
MM-RE-069




Fig. LLM
MM-RE-070
Bayonet Distal




Fig. LLM
MM-RE-071
Bayonet Mesial

Hylin



Fig. 1
MM-RE-072



Fig. 2
MM-RE-073




Fig. 3
MM-RE-074
Left




Fig. 4
MM-RE-075
Right



Bein



Fig. 1
MM-RE-076



Fig. 6
MM-RE-077



Fig. 7
MM-RE-078

Flohr




Fig. 1
MM-RE-079




Fig. 2
MM-RE-080
Left




Fig. 3
MM-RE-081
Right



Root Elevators

Schlemmer



Fig. 1
MM-RE-082
Left




Fig. 2
MM-RE-083
Right

Coleman / London Hospital



Fig. 20
MM-RE-086
Coleman



Fig. 21
MM-RE-087
London Hospital
Left




Fig. 22
MM-RE-088
London Hospital
Right

Schlemmer



Fig. 1
MM-RE-084
Left



Fig. 2
MM-RE-085
Right

Warwick-James



MM-RE-089
Straight



MM-RE-090
Left



MM-RE-091
Right

Root Tip Picks / Heidebrink



Fig. 2
MM-RE-093
Left




Fig. 1
MM-RE-092



Fig. 3
MM-RE-094
Right




Fig. 13/14
MM-RE-095
Double ended



Fig. 2
MM-RE-097
Left




Fig. 1
MM-RE-096






Fig. 3
MM-RE-098
Right



MM-RE-099
Lecluse Fig. 1



Potts




Fig. 6
MM-RE-100
Left



Fig. 7
MM-RE-101
Right

Winter




Fig. 1L
MM-RE-102
Left



Fig. 1R
MM-RE-103
Right



Fig. 11L
MM-RE-104
Left



Fig. 11R
MM-RE-105
Right




Fig. 12L
MM-RE-106
Left




Fig. 12R
MM-RE-107
Right



Fig. 13L
MM-RE-108
Left



Fig. 13R
MM-RE-109
Right




Fig. 13L
MM-RE-110
Left



Fig. 13R
MM-RE-111
Right

Walter F. Barry



Fig. 320
MM-RE-112
Left




Fig. 321
MM-RE-113
Right




Fig. 322
MM-RE-114
Left




Fig. 323
MM-RE-115
Right



Fig. 350
MM-RE-116
Left



Fig. 351
MM-RE-117
Right



Root Elevators



Fig. 1
MM-RE-001



Fig. 2
MM-RE-002



Fig. 3
MM-RE-003




Fig. 44
MM-RE-029
Left



Fig. 45
MM-RE-030
Right



Fig. 2
MM-RE-059
Left



Fig. 3
MM-RE-060
Right



MM-RE-118
Set of 7



Fig. 1
MM-RE-001



Fig. 2
MM-RE-002



Fig. 3
MM-RE-003




Fig. 44
MM-RE-029
Left



Fig. 45
MM-RE-030
Right



Fig. 2
MM-RE-059
Left



Fig. 3
MM-RE-060
Right



MM-RE-119
Set of 7



Fig. 1
MM-RE-001



Fig. 2
MM-RE-002



Fig. 3
MM-RE-003




Fig. 44
MM-RE-029
Left



Fig. 45
MM-RE-030
Right



Fig. 2
MM-RE-059
Left



Fig. 3
MM-RE-060
Right



MM-RE-120
Set of 7

Chompret



Fig. 1
MM-RE-121



Fig. 2
MM-RE-122



Fig. 3
MM-RE-123



Fig. 4
MM-RE-124





Root Elevators

Serrated




Fig. 1
MM-RE-125



Fig. 2
MM-RE-126



Fig. 1
MM-RE-127




Fig. 3
MM-RE-128



Fig. 4
MM-RE-129



Fig. 6
MM-RE-130



Fig. 7
MM-RE-131

Serrated Elevators MM-RE-125-131 Set of 7

Luxators



Fig. 1
MM-RE-132



Fig. 2
MM-RE-133




Fig. 3
MM-RE-134



Fig. 4
MM-RE-135




Fig. 5
MM-RE-136



Fig. 6
MM-RE-137




Fig. 7
MM-RE-138

Luxators MM-RE-132-138 Set of 7






Fig. 2
MM-RE-139
Flohr




Fig. 3
MM-RE-140
Flohr



Fig. 1
MM-RE-141
Bein



Fig. 6
MM-RE-142
Bein




Fig. 7
MM-RE-143
Bein

MM-RE-144 Set of 5

British Models




Fig. 1
MM-RE-145
Coupland




Fig. 2
MM-RE-146
Coupland






Fig. 3
MM-RE-147
Coupland



MM-RE-148
Warwick-James



MM-RE-149
Warwick-James



MM-RE-150
Warwick-James




Fig. 30
MM-RE-151
Cryer




Fig. 31
MM-RE-152
Cryer

MM-RE-153 Set of 8

Root Tip Elevators




Fig. 1
MM-RE-154
Heidebrink




Fig. 2
MM-RE-155
Heidebrink




Fig. 3
MM-RE-156
Heidebrink



Fig. 11
MM-RE-157
Bernhard




Fig. 6
MM-RE-158
Bernhard



Fig. 7
MM-RE-159
Bernhard

MM-RE-160 Set of 6

Page 84–110


- MOUTH GAGS
- TONGUE DEPRESSORS
- LIPS & CHEEK RETRACTORS
- RETRACTORS & FORCEPS
- BONE INSTRUMENTS
- IMPLANTOLOGY INSTRUMENTS
- PERIOSTEAL ELEVATORS
- PERIODONTAL INSTRUMENTS
- BONE CURETTES & FILES
- SURGERY INSTRUMENTS



Molt
MM-MG-001 10 cm
MM-MG-002 11 cm
MM-MG-003 14 cm



Roser-König
MM-MG-004 16 cm
MM-MG-005 19 cm



MM-MG-006
 Denhardt
 13 cm



MM-MG-007
 Fergusson - Ackland
 14.5 cm



MM-MG-008
 Heister
 Maunders Wood
 9.5 cm




MM-MG-009
 Obwegeser 11 cm
 Spreading 47 mm



MM-MG-010
 Heister
 13 cm

Mouth Gags



Davis-Boyle Tongue Depressors
With fixed suction tube

Fig. 1
MM-MG-011
110 x 28 mm



Fig. 2
MM-MG-012
95 x 27 mm



Fig. 3
MM-MG-013
80 x 27 mm



Fig. 4
MM-MG-014
65 x 26 mm



Fig. 5
MM-MG-015
30 x 20 mm



MM-MG-016
Davis-Boyle
Complete, consisting of frame
& 5 tongue depressors No. 1 to 5



MM-MG-017
Mc-Ivor
Frame for Mouth Gag



Mc-Ivor
Tongue Depressors

Fig. 1
MM-MG-018
100 x 25 mm

Fig. 2
MM-MG-019
85 x 25 mm

Fig. 3
MM-MG-020
75 x 25 mm

MM-MG-021
Mc-Ivor
Complete, with 3 tongue depressors



MM-MG-022
Dingmann
Complete, with frame, 3 tongue depressors,
& 2 adjustable lateral blades



Dingmann
Tongue Depressors
13.5 cm

Fig. 1
MM-MG-023
62 x 26 mm

Fig. 2
MM-MG-024
65 x 30 mm

Fig. 3
MM-MG-025
78 x 32 mm

Mouth Gags




Jennings
MM-MG-026 10.5 cm
MM-MG-027 12.5 cm
MM-MG-028 14.5 cm
MM-MG-029 16.5 cm



Whitehead Mouth Gag
MM-MG-030 11cm
MM-MG-031 13cm
MM-MG-032 15cm



Cheek Retractor
MM-MG-033 Small
MM-MG-034 Medium
MM-MG-035 Large



MM-TD-001
Hartmann
15 cm



MM-TD-002
Doyen Tongue Spatula
13.5 cm




MM-TD-003
14.5 cm




MM-TD-004
18 cm




MM-TD-005
Bruenings
19 cm




MM-TD-006
Tobold
15 cm




MM-TD-007
Tobold
16 cm



MM-TD-008
14 cm



MM-TD-009
15 cm



Wieder

Lips & Cheek Retractors



MM-RT-001
Obwegeser Chin
Retractor
16.5 cm



Obwegeser 16 cm
MM-RT-002 8 mm
MM-RT-003 10 mm



MM-RT-004 43 x 11 mm
MM-RT-005 57 x 11 mm
MM-RT-006 83 x 16 mm



Obwegeser
MM-RT-007
Spina-nasalis retractor
22 cm



Obwegeser
MM-RT-008
Ramus retractor for
separation of ramus & tissue
22 cm



MM-RT-009
Simplex
10.5 cm



30x19 mm

33x25 mm

MM-RT-010
Kilner
15 cm



MM-RT-011
Minnesota
14 cm



MM-RT-012
Cawood-Minnesota
15 cm




MM-RT-013
Austin



MM-RT-014
Hilger
Fits any simple stem mouth mirror handle



MM-RT-015
Hilger
with handle




MM-RT-016
20 cm

Finger protectors



Langenbeck

MM-RT-017 Small
MM-RT-018 Medium
MM-RT-019 Large



MM-RT-020
Oringer
11.5 cm



MM-RT-021
Obwegeser Cheek Retractor
12cm



MM-RT-022
Columbia
14.5 cm



Sternberg


MM-RT-023 12 cm
MM-RT-024 14 cm
MM-RT-025 15 cm



Vestibulum
MM-RT-026




MM-RT-027



MM-RT-028

Implants Retractors

Retractors & Forceps



Senn-Miller
16 cm



16 cm




16 cm




MM-RT-039
Collin Tongue holding forcep
17 cm



Young Tongue holding forcep
MM-RT-040 15 cm
MM-RT-041 17 cm




MM-RT-043
Farabeuf 12 cm
Pair



MM-RT-044
Parker-Langenbeck 21 cm
Pair




Kocher
21.5 cm



MM-RT-045
55 x 11 mm

MM-RT-046
40 x 11 mm



MM-RT-047
55 x 11 mm


MM-RT-048
40 x 11 mm




Middeldorph



MM-RT-049
14 x 17 mm
21.5 cm



MM-RT-050
20 x 22 mm
22 cm




MM-RT-051
26 x 30 mm
23.5 cm




■ Langenbeck
22 cm




● Langenbeck
22 cm




■ MM-RT-052
● MM-RT-053




■ MM-RT-054
● MM-RT-055



■ MM-RT-056
● MM-RT-057





■ MM-RT-058
● MM-RT-059




■ MM-RT-060
● MM-RT-061

Bone Instruments


Chisels & Gouges




Osteotomes



Pocket Markers



Elevators & Raspatories



PDL Elevators



MM-IM-001 Distal, 2.5 mm Micro Edge



MM-IM-002 Mesial, 2.5 mm Micro Edge



MM-IM-003 Distal, 4 mm Micro Edge



MM-IM-004 Mesial 4 mm Micro Edge



MM-IM-005 2.5 mm Micro Edge, Straight



MM-IM-006 2.5 mm Micro Edge, Curved




MM-IM-007 3.4 mm Straight

MM-IM-008 Set Of 7 Pieces



MM-IM-009 Distal, 2.5 mm Micro Edge



MM-IM-010 Distal, 2.5 mm Micro Edge



MM-IM-011 Mesial, 2.5 mm Micro Edge



MM-IM-012 Mesial 4 mm Micro Edge



MM-IM-013 2.5 mm Micro Edge, Straight





MM-IM-014 2.5 mm Micro Edge, Curved




MM-IM-015 3.4 mm Straight

MM-IM-016 Set Of 7 Pieces

Osteotomes





MM-IM-031
Bone Collector Chisel Straight



MM-IM-033
Straight blade for bone collector



MM-IM-032
Bone Collector Chisel Curved



MM-IM-034
Curved blade for bone collector



MM-IM-035
Ratchet Screw driver handle
Hexagonal head



MM-IM-036
Bone Morselizer with Titanium Grinding



MM-IM-037
Bone Mixing Basin
Dia 35 mm



MM-IM-038
Manual Easy Bone Mill



MM-IM-039
Bone Craker Titanium

Sinus Lift Currettes



Sinus Lift Elevators

Hollow handle



Hollow handle



TiN Coated



MM-IM-056
Fig. 1



MM-IM-064
Fig. 1



MM-IM-057
Fig. 2



MM-IM-065
Fig. 2



MM-IM-058
Fig. 3



MM-IM-066
Fig. 3



MM-IM-059
Fig. 4



MM-IM-067
Fig. 4



MM-IM-060
Fig. 5



MM-IM-068
Fig. 5



MM-IM-061
Fig. 6



MM-IM-069
Heights Fig. 6



MM-IM-062
Fig. 7



MM-IM-070
Fig. 7



MM-IM-063
Fig. 8



MM-IM-071
Fig. 8

Bone Collector




MM-IM-072

MM-IM-073

Disposable stainless steel cartridge
Pack of 10 pieces

Applicators for bone grafting material



TIN coated

MM-IM-074 Applicator



MM-IM-075 Applicator & compactor

Compactors



MM-IM-076 Hollow handle



MM-IM-077 Hollow handle

TIN coated



MM-IM-078 Hollow handle



MM-IM-079 Hollow handle

Membrane placement instruments




MM-IM-080




MM-IM-081
11 cm

Osteotomes




Dilator



Spreader

A-Form
Ø 2 mm Ø 3 mm Ø 3.3 mm Ø 4 mm




MM-IM-082
14.5 cm

MM-IM-083
14.5 cm

MM-IM-084
14.5 cm

MM-IM-085
14.5 cm

D-Form
Convex Tips
Ø 3 mm Ø 4 mm Ø 5 mm



MM-IM-086
14.5 cm

MM-IM-087
14.5 cm


MM-IM-088
14.5 cm



14.5 mm



FREER Osteotomos



FREER Osteotomos

Ø 1.6 mm Ø 1.9 mm Ø 2.8 mm Ø 3.1 mm Ø 3.9 mm



MM-IM-089

MM-IM-090

MM-IM-091

MM-IM-092

MM-IM-093

Ø 3.0 mm Ø 3.3 mm Ø 4.0 mm





MM-IM-094

MM-IM-095

MM-IM-096

Osteotomes






MM-IM-119
Trephine Drills Stand
Complete with drills

Tissue Punches



MM-IM-120
Trephine Drills Stand
Complete with drills



MM-IM-121
3 mm

MM-IM-124
4 mm

MM-IM-122
5 mm

MM-IM-123
8 mm

Keyes



MM-IM-125
3 mm

MM-IM-126
4 mm

MM-IM-127
5 mm

MM-IM-128
6 mm

Keyes

Bone Thickness Gauges



Fig. 1
Castroviejo
MM-IM-129



Fig. 2
Castroviejo
MM-IM-130




Fig. 3
Castroviejo
MM-IM-131



Fig. 4
Castroviejo
MM-IM-132



Fig. 5
Castroviejo
MM-IM-133



MM-IM-134
Bone Thickness Gauge



MM-IM-135
Bone Caliper
15 cm



MM-IM-136
Sinus Rongeur
Delicate bone cutting forceps
for opening the sinus window
4 mm bite, TiN coated

Bone Injectors



Fig. 1
MM-IM-137



Fig. 2
MM-IM-138



Fig. 3
MM-IM-139




Fig. 4
MM-IM-140



TiN coated bone block clamps

To hold & stabilize block graft against recipient site while placing screws.
Slotted tips act as a guide for precise drilling.



MM-IM-141
Posterior

MM-IM-142
Anterior

Depth measuring probes



MM-IM-143
Titanium tip, silicone handle & laser marked

Periotome



Fig. 1 MM-IM-144

Implant periosteal elevators & raspatories



Fig. 1 MM-IM-145



Fig. 2 MM-IM-146 Sharp-Blunt



Fig. 1 MM-IM-147 Atraumatic papilla elevator

Mallets



MM-IM-148
Partsch
Stainless Steel head
19 cm



MM-IM-149
Partsch
Lead filled head
19 cm



MM-IM-150
Mead 19 cm
MM-IM-151
Spare plastic disks, pair



MM-IM-152
18.5 cm
MM-IM-153
Spare plastic disks, pair



MM-IM-154
19 cm
MM-IM-155
Spare plastic disks, pair



MM-IM-156
20 cm
MM-IM-157
Spare plastic disks, pair

Periosteal Elevators



Fig. 15 **MM-PE-001** Freer



Fig. 16 **MM-PE-002** Freer




Fig. 14 **MM-PE-003** Goldman-Fox




Fig. 20 **MM-PE-004** Hirschfield



Fig. 1 **MM-PE-005** Hopkins




Fig. 2 **MM-PE-006** Hopkins




Fig. 3 **MM-PE-007** Mead/West



Fig. 9 **MM-PE-008** Molt




Fig. 3 **MM-PE-009** Prichard




Fig. 9 **MM-PE-010** Molt



Fig. 1 MM-PE-011 Busar Blue



Fig. 3 MM-PE-012 Prichard Blue

Implant periosteal elevators



Fig. 1 MM-PE-013 Busar



Fig. 9A MM-PE-014 Allen



Fig. 25 MM-PE-015 17 cm

Gingivectomy Knives

Hollow handle




Fig. 15/16 MM-PE-016 Kirkland



Fig. 1/2 MM-PE-017 Orban

Solid handle



Fig. 15/16 MM-PE-018 Kirkland




Fig. 1/2 MM-PE-019 Orban








Fig. 1 MM-PE-020 Orban





Fig. 2 MM-PE-021 Orban





Bone Curettes





	Solid handle	Hollow handle
Fig. 1	MM-BC-014	MM-BC-018
Fig. 2	MM-BC-015	MM-BC-019
Fig. 3	MM-BC-016	MM-BC-020
Fig. 4	MM-BC-017	MM-BC-021



	Solid handle	Hollow handle
Fig. 85	MM-BC-022	MM-BC-026
Fig. 86	MM-BC-023	MM-BC-027
Fig. 87	MM-BC-024	MM-BC-028
Fig. 88	MM-BC-025	MM-BC-029

Bone Curettes & Files

Bone curettes



Molt
16.5 cm




Fig. 1
MM-BC-030




Fig. 2
MM-BC-031




Fig. 4
MM-BC-032





Fig. 5
MM-BC-033



Fig. 6
MM-BC-034



Williger
14 cm




Fig. 00/0
MM-BC-035




Fig. 0/1
MM-BC-036




Fig. 1/2
MM-BC-037

Bone Files



Fig. 10
MM-BF-001



Fig. 12A
MM-BF-002



Fig. 33
MM-BF-003



Fig. 12
MM-BF-004
Howard




Fig. 45
MM-BF-005
Miller




Fig. 64
MM-BF-006
Miller




Skin Hooks



Probes




Self-Retaining Retractors



Suture Forceps



Suction Tubes





Page 111–115

- SELF & MANUAL ASPIRATING SYRINGES
- INTRALIGAMENTAL SYRINGES



MM-SY-001
Type - A
Metric Thread
1.8 ml



MM-SY-002
Type - C
Metric Thread
1.8 ml


Manual Aspirating Syringes



MM-SY-003
Type - A
Metric Thread
1.8 ml



MM-SY-004
Type - C
Metric Thread
1.8 ml



MM-SY-005
P - Ject
Metric Thread
1.8 ml



MM-SY-006
Type - A
Metric Thread
1.8 ml




MM-SY-007
Euro
Metric Thread
1.8 ml


Manual Aspirating Syringes




MM-SY-008
Euro
Metric Thread
1.8 ml




MM-SY-009
T - Handle
Metric Thread
1.8 ml



MM-SY-010
3 - Ring
Metric Thread
1.8 ml



MM-SY-011
Euro breech loading, Ring handle
Metric Thread
1.8 ml



MM-SY-012
Euro breech loading, T - handle
Metric Thread
1.8 ml




MM-SY-013
Euro
Metric Thread
1.8 ml




MM-SY-014
3 - Ring
Metric Thread
1.8 ml



MM-SY-015
2 - Ring
Metric Thread
1.8 ml




MM-SY-016
Euro breech loading, Ring handle
Metric Thread
1.8 ml



MM-SY-017
Euro breech loading, T - handle
Metric Thread
1.8 ml

Intraligamental Syringes



MM-SY-018
Curved
Metric Thread
1.8 ml



MM-SY-019
Curved
Metric Thread
1.8 ml



MM-SY-020
Straight
Metric Thread
1.8 ml



MM-SY-021
Swedish Type
Metric Thread
1.8 ml



Page

116–129

- RUBBER DAM INSTRUMENTS
- MATRIX RETAINERS & RUBBER DAM FRAMES
- LIGATURE & SEPARATORS PLACING PLIERS
- ELASTIC LIGATURE PLACING FORCEPS
- ELASTIC APPLICATORS & LIGATURE TUCKERS
- BRACKET POSITIONING INSTRUMENTS
- BRACKET & MEASURING INSTRUMENTS
- CROWN REMOVERS & FORCEPS FOR ARTICULATING PAPER
- CROWN SPREADERS
- CROWN INSTRUMENTS
- LABORATORY INSTRUMENTS



MM-RD-001
Brewer



MM-RD-002
Ivory



MM-RD-003
Ivory



MM-RD-004
Matrix Plier
17 cm



MM-RD-005
Stokes



MM-RD-006
Stokes / Palmer



MM-RD-007
Ivory

Rubber Dam Punches




Fig.1
MM-RD-008
Ainsworth



Fig.2
MM-RD-009
Ainsworth


Rubber Dam Instruments

RIGHT TEETH


MM-RD-010 211 Winged	MM-RD-011 2A Winged	MM-RD-012 212 Wingless	MM-RD-013 29 Wingless
MM-RD-020 4 Winged	MM-RD-021 205 Winged	MM-RD-022 5 Winged	MM-RD-023 8A Winged
MM-RD-024 201 Winged	MM-RD-025 14 Winged	MM-RD-026 51 Wingless	MM-RD-027 26 Wingless
MM-RD-028 W8A Wingless	MM-RD-029 18 Wingless		
MM-RD-030 1 Winged	MM-RD-031 1A Winged	MM-RD-032 207 Winged	MM-RD-033 208 Winged
MM-RD-034 2A Winged	MM-RD-035 22 Wingless	MM-RD-036 27 Wingless	MM-RD-037 W2A Wingless
MM-RD-038 11 Winged	MM-RD-039 10 Winged	MM-RD-040 3 Winged	MM-RD-041 203 Winged
MM-RD-042 204 Winged	MM-RD-043 202 Winged	MM-RD-044 201 Winged	
MM-RD-045 211 Winged	MM-RD-046 2 Winged	MM-RD-047 22 Wingless	

LEFT TEETH

MM-RD-055 208 Winged	MM-RD-056 207 Winged	MM-RD-057 206 Winged	MM-RD-058 22 Wingless	MM-RD-059 27 Wingless	MM-RD-060 29 Wingless
MM-RD-061 W3 Wingless	MM-RD-062 31 Wingless	MM-RD-063 30 Wingless	MM-RD-064 28 Wingless	MM-RD-065 25 Wingless	MM-RD-066 15 Wingless
MM-RD-067 205 Winged	MM-RD-068 5 Winged	MM-RD-069 8A Winged	MM-RD-070 201 Winged	MM-RD-071 14A Winged	MM-RD-072 4 Wingless
MM-RD-073 8 Wingless	MM-RD-074 28 Wingless	MM-RD-075 24 Wingless			
MM-RD-076 207 Winged	MM-RD-077 209 Winged	MM-RD-078 0 Winged	MM-RD-079 00 Winged	MM-RD-080 2A Winged	MM-RD-081 1 Winged
MM-RD-082 1A Winged					
MM-RD-083 200 Winged	MM-RD-084 202 Winged	MM-RD-085 14 Winged	MM-RD-086 56 Winged	MM-RD-087 7 Winged	MM-RD-088 7A Winged
MM-RD-089 14A Winged	MM-RD-090 56A Winged				
MM-RD-091 29 Wingless	MM-RD-092 22 Wingless	MM-RD-093 27 Wingless	MM-RD-094 W2A Wingless		




UPPER
PERMANENT
LOWER




MM-RD-095
Rubber Dam Clamps Tray


Rubber Dam Frames



MM-RD-096 13 cm
MM-RD-097 15 cm
Young



MM-RD-098
15 cm



MM-RD-099

Matrix Retainers



MM-MR-001
with cheek protector
5 mm



Fig. 8
MM-MR-002
Ivory
5 mm



Fig. 8
MM-MR-003
Ivory
6 mm



Fig. 8
MM-MR-004
Ivory
7 mm



MM-MR-005
Tofflemire junior



MM-MR-006
Tofflemire senior



MM-MR-007
Tofflemire universal



MM-MR-008
Nystrom
Left



MM-MR-009
Nystrom
Right



Fig. 12A
MM-MR-010
Ivory









Fig. 12
MM-MR-011
Ivory



MM-MR-012
Strip clamp
straight



Ligature & Separators Placing Pliers






MM-EL-007
Str 12 cm



MM-EL-008
Cvd 12 cm



MM-EL-009
14 cm



MM-EL-010
12 cm



MM-EL-011
14 cm



MM-EL-012
Elastomeric ligature
removing tweezer



Ligature plier
14 cm



Fig.1
MM-EL-013



Fig.2
MM-EL-014



Fig.3
MM-EL-015



Smaha
Ligature plier
14 cm



MM-EL-016
Str



MM-EL-017
Cvd



MM-EL-018
Ligature plier
14 cm



With tungsten
carbide insert

Elastic Applicators & Ligature Tuckers

Elastic Ligature Applicators

For easy placement of elastic ligatures under bracket tie wings



MM-EL-019



MM-EL-020

Elastomeric Ligature Removing Instruments

Multi-purpose instruments, ligature tucker & elastomeric ligature remover



MM-EL-021



MM-EL-022

Ligature Tuckers



MM-EL-023



MM-EL-024



MM-EL-025

Ligature Tuckers

Solid handle



MM-EL-026



MM-EL-027



MM-EL-028

Distal End Bending Instruments



MM-BP-001
For wire upto .22" (0.56 mm)

MM-BP-002
For wire upto .30" (0.76 mm)



MM-BP-003
For wire upto .22" (0.56 mm) Green

MM-BP-004
For wire upto .30" (0.76 mm) Green



Lingual Torquing Key

Orange **MM-BP-005**



NiTi Ice Instrument 11.5 cm

MM-BP-006
Turns NiTi wire soft and flexible immediately for faster wire placement



Bracket Placer for Posterior Teeth



MM-BP-009

Bracket Positioning Gauges
Anodized Aluminium



MM-BP-010 For .018" slot brackets



MM-BP-011 For .022" slot brackets



MM-BP-012
Boone
Universal

Adhesive Removers



MM-BP-007 **MM-BP-008**
Yellow

Bracket Placing Tweezers



MM-BP-013 13.5 cm



MM-BP-014 13.5 cm



MM-BP-015 11 cm

Bracket Tweezers with Slot Aligner



MM-BP-016 13.5 cm

Fig. 1



MM-BP-017 13.5 cm

Fig. 2



MM-BP-018 13.5 cm

Fig. 3

Bracket Measuring Instruments



MM-BP-019
Band Seater
15 cm




MM-BP-020
Mershon band pusher
14 cm



MM-BP-021
Band seating

Scaler & band pushing Instruments



MM-BP-022



MM-BP-023
Brown

Reverse action buccal tube tweezer



MM-BP-024 12.5 cm

MM-BP-025 13.5 cm

Buccal tube tweezer



MM-BP-026 12 cm

MM-BP-027 13.5 cm
Keat

Measuring Instruments




Bernklau
MM-BP-028 Orthodontic ruler
Perfect for all diagnostic measurements,
both intraoral and on the model



Ruler stainless steel
MM-BP-029 15 cm / 6"
MM-BP-030 30 cm / 12"




Beerendonk
MM-BP-031
Dental vernier with long & fine tips for intraoral measurements



MM-BP-032
Double sided, for inner & outer measurements



MM-BP-036
Movable Heads Bracket Positioning Gauge



Orthodontic bow compass
MM-BP-033

Turrets


For forming archwires, without torque,
for wires .016" to .022" (0.41mm to 0.56mm)




MM-BP-034
Made of Stainless steel



MM-BP-035
Made of anodized Aluminum
Disinfectable




MM-CR-001
Contents: 1 Crown Remover, Mounting Kit (Wrench & Key), 3 Tips




MM-CR-006
Contents: 1 Crown Remover, Mounting Kit (Wrench & Key), 3 Tips

Morrel Manual Crown Remover




MM-CR-011
Complete with 3 Tips

Morrel Manual Crown Remover with Spring



MM-CR-015
Complete with 3 Tips

Semi-Automatic Crown & Bridge Remover



MM-CR-019
Contents: 1 Crown Remover, Tightening Key, 3 Tips

Forceps For Articulating Paper



MM-CR-024
Miller
15 cm

MM-CR-025
Miller
15 cm

MM-CR-026
15 cm

Crown Removers & Crown Spreaders

manual crown remover double ended



Fig. 1
MM-CR-028



Fig. 2
MM-CR-029





Fig. 3
MM-CR-030



MM-CR-027
Complete with 3 tips



MM-CR-031
Wynman
Crown Plier
with tungsten carbide inserts



MM-CR-032
Furrer Crown Plier
16 cm



MM-CR-033
Crown spreader



Fig. 3
MM-CR-034
Crown Spreader



Crown Spreader






Fig. 0
MM-CR-035



Fig. 1
MM-CR-036



Fig. 2
MM-CR-037



Furnace Tongs



MM-LI-001
31 cm



MM-LI-002
40 cm



MM-LI-003
56 cm



MM-LI-004 20 cm Str
MM-LI-005 20 cm Cvd
Plaster Cutter

Calipers



Fig.1
MM-LI-006
Boley
Measuring range
0-100 mm




Fig.2
MM-LI-007
Boley
Measuring range
0-100 mm



MM-LI-008
Iwanson for Wax
Measuring range
0-10 mm



MM-LI-009
Iwanson for Metal
Measuring range
0-10 mm




MM-LI-010
Micro Boley Gauge

Modeling saw & blades



MM-LI-011

Fret saw & blades



MM-LI-017

Saw blades for metal Pack of 12 pcs



MM-LI-012

0.25 x 3 mm Pack of 12 pcs



MM-LI-013

0.25 x 5 mm Pack of 12 pcs



MM-LI-014

0.12 x 5 mm Pack of 12 pcs



MM-LI-015

Soldering Tweezer
Str 16 cm



MM-LI-016

Soldering Tweezer
Cvd 16 cm



Fig.1

MM-LI-019

Riveting Mallet
22 cm



Fig.2

MM-LI-020

Riveting Mallet
22 cm



MM-LI-021

Anvil
Carbon Steel

Bib chains



Fig. 1
MM-LI-022
Ball chain



MM-LI-023
Plastic chain

X-Ray holders



Fig. 1
MM-LI-024



Fig. 6
MM-LI-025
6 Clips



Fig. 10
MM-LI-026
10 Clips




Fig. 14
MM-LI-027
14 Clips

Page
130–131

○ ARTICULATORS

Articulators



MM-AR-001
Articulator ASA W/out Bearing Rod
Anatomical Articulator With Incisive
Pin, 4 Mounting Plates And Supporting
Back-pin. Available In Assorted Colors




MM-AR-002
ASA Articulator With Lateral Moment.
Made Of Aluminum Die Casting With
2 Mounting Plates. Also Available In
Vibratory Polished Process & Assorted colors



MM-AR-003
Keystone Deluxe Magnetic Denture
Articulator With Plates, Interchangeable,
Incisal Pin Included. Magnetic Split Casting
Ensures Easy & Exact Setting, Lateral
Excursions



MM-AR-004
Keystone Deluxe Universal Type-A Articulator.
Interchangeable, Included Incisal Pin, Lateral
Excursions



MM-AR-005
Anatomical Articulator Universal, Realized
For Easy Visibility To The Lingual Surface
& For A Perfect Stability Also In Reserved
Position. Also available in Assorted Colors



MM-AR-006
Chromium Plated Fork
Fork Crowned Bridge Nickel
Chromium Plated Articulator



MM-AR-007
Full Arch Hinged Articulator
Chromium Plated Hinge Articulator
For Mounting Full Arch Cases. Low
Back. Gypsum-free Articulation Enables
Clean & Accurate Mounting & Dismounting
In 30 Seconds



MM-AR-008
Ray Foster 57 Chrome Articulator
Chrome Brass Hinge Crowned Bridge
Articulator For Mounting Quadrant Cases.
Crown & Bridge Articulator Is Popular Ray
Foster # 57 Design, Eliminates The Need
To Break Off Models



MM-AR-009
Occlusor, Wire Articulator, Available
In 2 Sizes, Small & Large



MM-AR-010
Professional Occlusor, Vibratory
Polished Articulator.



MM-AR-011
Dental Articulator 2x-cresant Brand
Full Denture Chromium Plated Articulator.
Provides Additional Vertical Dimension For
Acceptance Of Higher Models. 49 Wide At
The Joint & Has A Vertical Span Between
Parallel Bows Of 25/8"



MM-AR-012
Hinge Articulator Full Arch
Plaine Line Nickle Chromium
Plated Articulator



MM-AR-013
Magic Art-I Articulator Lightweight



MM-AR-014
Mortimer-I Oculader Articulator




MM-AR-015
Mortimer-II Oculader Articulator

Page

132–149

- ORTHODONTIC PLIERS
- WIRE CUTTERS WITH TUNGSTEN CARBIDE INSERTS
- BONE CUTTERS
- INSTRUMENTS FOR "FIXED APPLIANCE TECHNIQUE"
- BONE RONGEURS



MM-OP-001
Flat nose wire bending plier
13 cm



MM-OP-002
Peeso
Round nose wire bending plier
13 cm



MM-OP-003
Peeso
Crown & collar cotouring plier
13 cm



MM-OP-004
Aderer-Mini
11.5 cm




MM-OP-005
Aderer-Maxi
12 cm



MM-OP-006
Schwarz graduated
Optical style
13 cm



MM-OP-007
Schwarz
Optical style, round &
concave beaks 14 cm










MM-OP-008
Schwarz
Optical style, round &
concave beaks 14 cm







MM-OP-009
Mckellcps
Optical style, round &
concave beaks 16 cm

Orthodontic pliers





Replacement pads
MM-OP-023 4.8mm (1/4"), 2 pieces
MM-OP-024 6.3mm (3/16"), 2 pieces



Orthodontic pliers



MM-OP-031
Replacement tip stainless
440 fcr Omega or Tweed
loop forming pliers.
Tip diameters: 1.14mm, 1.52mm
& 1.91 mm




MM-OP-032
Omega 13 cm
Loop forming pliers with replaceable
stainless steel tip, three step tip has
1.14, 1.52 & 1.91 mm diameters



MM-OP-033
Tweed
Loop forming plier
13 cm



MM-OP-034
Angled




How
standard
4 mm wide tips
13 cm



MM-OP-035
Straight



MM-OP-036
Band Seating
Plier 13 cm



MM-OP-037
Nance
Loop forming plier
with four steps beak
13 cm



MM-OP-038
Tweed
Loop forming plier
13 cm









MM-OP-039
Weingart Mini
13 cm



MM-OP-040
Weingart
13 cm



MM-OP-041
Weingart Narrow
Narrow long tips for better
access to reach hard areas
13 cm



Orthodontic pliers



MM-OP-051
Adams Mini
12 cm



MM-OP-052
Adams Maxi
12.5 cm



MM-OP-053
Round nose plier
Fine pattern
13 cm



MM-OP-054
Adams Mini
13 cm



MM-OP-055
Adams Maxi
13 cm



MM-OP-056
Light wire bending
& loop forming Plier
13 cm



MM-OP-057
Angel
Bird beak plier
13 cm



MM-OP-058
Angel
Round nose wire
bending plier
13 cm



MM-OP-059
Angel
Round nose wire bending
& loop forming plier
13 cm







MM-OP-060
Serrated Laboratory
plier 14 cm



MM-OP-061
Smooth Laboratory
plier 14 cm



MM-OP-062
Laboratory plier
14.5 cm




MM-OP-063
Ribbon Arch plier
13 cm




MM-OP-064
Cone & Socket plier
13 cm



MM-OP-065
Goslee contouring
plier 13 cm



MM-OP-066
NiTi Cinch back
plier



MM-OP-067
Triplex Contouring
plier
14.5 cm



MM-OP-068
Wire bending
plier
14 cm

Orthodontic pliers



MM-OP-069
De La Rosa
Graduated
13 cm



MM-OP-070
De La Rosa
Smooth
13 cm



MM-OP-071
Wire bending plier
13 cm



MM-OP-072
Andresen wire
bending plier
13 cm




MM-OP-073
Hawley retainer
Plier
13 cm




MM-OP-074

Universal-Dolphin style
Wire bending & cutting plier
16 cm

● Bends: 0.9 mm (.036")
Cuts: 0.7 mm (.028")


MM-OP-078

Nance
Adams clasp forming plier
13 cm

Max. 0.70 mm (.028")


MM-OP-075

Universal-Taper style
Wire bending & cutting plier
16 cm

● Bends: 0.9 mm (.036")
Cuts: 0.7 mm (.028")


MM-OP-079

Schwarz
Arrow clasp forming plier
13 cm

Max. 0.70 mm (.028")


MM-OP-076

Waldasachs
Wire bending & cutting plier
15 cm

● Bends: 0.9 mm (.036")
Cuts: 0.7 mm (.028")


MM-OP-080

Schwarz
Arrow clasp forming plier
13 cm

Max. 0.70 mm (.028")


MM-OP-077

Universal
parallel jaw Plier
with side cutter

Max. 1.80 mm (.071")


MM-OP-081

Magnum
Universal facebow loop
forming & cutting plier

Max. 1.80 mm (.071")

Wire Cutters with Tungsten Carbide Inserts



MM-OP-082
Distal end cutter
with long handle
14 cm

- Max. 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" x .025")
- ⊗ Max. 0.36 mm to 0.53 mm (.0155" to .0215")
- ⊗ Max. 0.53 mm x 0.64 mm (.021" x .025")




MM-OP-083
Distal end cutter
with long handle
14 cm

- Max. 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" x .025")
- ⊗ Max. 0.36 mm to 0.53 mm (.0155" to .0215")
- ⊗ Max. 0.53 mm x 0.64 mm (.021" x .025")



MM-OP-084
Distal end cutter
14 cm

- Max. 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" x .025")
- ⊗ Max. 0.36 mm to 0.53 mm (.0155" to .0215")
- ⊗ Max. 0.53 mm x 0.64 mm (.021" x .025")



MM-OP-085
soft wire & ligature cutter
13 cm

.25 mm to .38 mm (.010" to .015")



MM-OP-086
soft wire & ligature cutter
13 cm

.25 mm to .38 mm (.010" to .015")



MM-OP-087
soft wire & ligature cutter
45° angled
13 cm

.25 mm to .38 mm (.010" to .015")

Heavy Box Joint Construction Cutters



MM-OP-088
Hard wire cutter, straight
14 cm

- Max. 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.53 mm x 0.71 mm (.018" x .025")



MM-OP-089
Hard wire cutter, 15° angled
14 cm

- Max. 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.53 mm x 0.71 mm (.022" x .028")

With Tungsten Carbide Inserts



MM-OP-090
How Standard
4 mm wide tips
13 cm



MM-OP-091
Angled

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



MM-OP-092
Straight

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



MM-OP-093
Weingart Mini
13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



MM-OP-094
Weingart
13 cm

● Max. 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



MM-OP-095
Soft wire & ligature cutter
13 cm



□ .25 mm to .38 mm (.010" to .015")



MM-OP-096
Goldman-Fox
10 cm



Liston Bone cutter

MM-OP-097 14 cm Str **MM-OP-101** 14 cm Cvd
MM-OP-098 17 cm Str **MM-OP-102** 17 cm Cvd
MM-OP-099 19 cm Str **MM-OP-103** 19 cm Cvd
MM-OP-100 22 cm Str **MM-OP-104** 22 cm Cvd

Instruments for "Fixed Appliance Technique"

Standard distal end cutter with safety hold

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



MM-OP-105

13 cm

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.56 mm x 0.71 mm (.022" x .028")

Soft wire & ligature cutter, Mini

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



MM-OP-109

13 cm

- └ 25 mm to .38 mm (.010" to .015")

Mini head distal end cutter with safety hold

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



MM-OP-106

13 cm

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" x .025")
- ⊗ 0.38 mm to 0.53 mm (.0155" to .0215")
- ⊗ Max. 0.53 mm x 0.64 mm (.021" x .025")

Soft wire and ligature cutter, standard

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



MM-OP-110

13 cm

- └ 25 mm to .38 mm (.010" to .015")

Distal end cutter with safety hold, mini head, long handles

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



MM-OP-107

13 cm

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" x .025")
- ⊗ 0.38 mm to 0.53 mm (.0155" to .0215")
- ⊗ Max. 0.53 mm x 0.64 mm (.021" x .025")

Soft wire & ligature cutter, 15° angled

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



MM-OP-111

13 cm

- └ 25 mm to .38 mm (.010" to .015")

Flush-cut distal end cutter with safety hold

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



MM-OP-108

13 cm

- 0.30 mm to 0.50 mm (.012" to .020")
- Max. 0.46 mm x 0.64 mm (.018" x .025")

Soft wire & ligature cutter, 45° angled

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



MM-OP-112

13 cm

- └ 25 mm to .38 mm (.010" to .015")

Utility pliers



MM-OP-113
Weingart-Mini
13 cm
● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



MM-OP-114
Weingart
13 cm
● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Utility pliers



MM-OP-115
How
Straight 13 cm
● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")



MM-OP-116
How
Angled 13 cm
● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Wire bending & loop forming pliers



MM-OP-117
Angel (Bird beak)
13 cm
● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Light wire bending & loop forming pliers



MM-OP-118
13 cm
● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.41 mm x 0.56 mm (.016" x .022")

Light wire bending pliers with cutter Fitted with hard metal tips



MM-OP-119
13 cm
● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.46 mm x 0.56 mm (.018" x .022")

Light wire bending & loop forming pliers with cutter



MM-OP-120
13 cm
● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.46 mm x 0.56 mm (.018" x .022")

Light wire bending & loop forming pliers, grooved tip



MM-OP-121
Tweed
13 cm
● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.46 mm x 0.56 mm (.018" x .022")

Instruments for "Fixed Appliance Technique"

Hard wire cutter, straight

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



MM-OP-122
13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Bracket removing pliers angled



MM-OP-126
13 cm

Hard wire cutter, 15° angled

Non-autoclavable



MM-OP-123
13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Bracket removing pliers straight



MM-OP-127
13 cm

Hard wire cutter, 15° angled

Fitted with un-breakable, corrosion resistant & durable stainless steel tips



MM-OP-124
13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Bracket removing pliers with plastic tip



MM-OP-128
13 cm

Band removing pliers



MM-OP-125
13 cm

Adhesive removing pliers



MM-OP-129
13 cm

MM-OP-130 Replacement metal tip

Loop forming pliers with cutter



MM-OP-131

Tweed
13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.53 mm x 0.64 mm (.021" x .025")

Archwire forming pliers, concave & convex tips



MM-OP-135

Hollow-Chop
13 cm

● 0.30 mm to 0.71 mm (.012" to .028")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Light wire bending pliers



MM-OP-132

Jarabak
13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.46 mm x 0.56 mm (.018" x .022")

Multi-purpose pliers (bird beak, De la rosa & cutter in one)



MM-OP-136

13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.46 mm x 0.56 mm (.018" x .022")

3-jaw wire bending pliers



MM-OP-133

Aderer
13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Distal end bending pliers




MM-OP-137

(Cinch back / tie back)
13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Ligature forming pliers



MM-OP-134

13 cm

☞ .25 mm to .38 mm (.010" to .015")

Distal end bending pliers



MM-OP-138

(Hammer head)
13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Instruments for "Fixed Appliance Technique"

Loop forming pliers



MM-OP-139

Nance
13 cm

● 0.36 mm to 0.71 mm (.014" to .028")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Torquing pliers set



0.18" 0.22"

MM-OP-143

13 cm

Detailing step pliers / Bayonet bending pliers - 1mm



MM-OP-140

13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Rectangular arch forming & torquing pliers



MM-OP-144

Tweed / angel
13 cm

● 0.30 mm to 0.71 mm (.012" to .028")
■ Max. 0.56 mm x 0.71 mm (.022" x .028")

Stop pliers V- Bend - 1mm



MM-OP-141

13 cm

● 0.30 mm to 0.50 mm (.012" to .020")
■ Max. 0.43 mm x 0.56 mm (.017" x .022")

Surgical ball hook crimping pliers, angled



MM-OP-145

13 cm

Lingual / palatal arch forming pliers



MM-OP-142

13 cm









● 0.76 mm (.030") / 0.91 mm (.036")

Surgical ball hook crimping pliers, straight




MM-OP-146

13 cm




BoneRongeurs



MM-RB-011

MM-RB-012

Friedmann-Mini
12.5 cm



MM-RB-013
Friedmann
14 cm



MM-RB-014



MM-RB-015




MM-RB-016




Friedmann-Micro
15 cm




Blumenthal
14.5 cm



MM-RB-017



MM-RB-018



MM-RB-019



MM-RB-020
Bayer
18 cm


Page 150–154

○ IMPRESSION TRAYS

Perforated without Rim




Lower



No.	Size	A mm	B mm	Code
L1	XXS	55	47	MM-IT-001
L2	XS	60	49	MM-IT-002
L3	S	72	54	MM-IT-003
L4	M	74	57	MM-IT-004
L5	L	76	59	MM-IT-005
L6	XL	78	61	MM-IT-006



Upper




No.	Size	A mm	B mm	Code
U1	XXS	61	47	MM-IT-007
U2	XS	64	51	MM-IT-008
U3	S	66	53	MM-IT-009
U4	M	69	57	MM-IT-010
U5	L	71	60	MM-IT-011
U6	XL	73	62	MM-IT-012

Set	Contents	Code
6 Pcs	L3,L4,L5,U3,U4,U5	MM-IT-013
8 Pcs	L3,L4,L5,L6,U3,U4,U5,U6	MM-IT-014
12 Pcs	L1,L2,L3,L4,L5,L6 U1,U2,U3,U4,U5,U6	MM-IT-015

Perforated with Rim




Lower



No.	Size	A mm	B mm	Code
L1	XXS	55	47	MM-IT-016
L2	XS	60	49	MM-IT-017
L3	S	72	54	MM-IT-018
L4	M	74	57	MM-IT-019
L5	L	76	59	MM-IT-020
L6	XL	78	61	MM-IT-021



Upper




No.	Size	A mm	B mm	Code
U1	XXS	61	47	MM-IT-022
U2	XS	64	51	MM-IT-023
U3	S	66	53	MM-IT-024
U4	M	69	57	MM-IT-025
U5	L	71	60	MM-IT-026
U6	XL	73	62	MM-IT-027

Set	Contents	Code
6 Pcs	L3,L4,L5,U3,U4,U5	MM-IT-028
8 Pcs	L3,L4,L5,L6,U3,U4,U5,U6	MM-IT-029
12 Pcs	L1,L2,L3,L4,L5,L6 U1,U2,U3,U4,U5,U6	MM-IT-030

Solid without Rim




Lower



No.	Size	A mm	B mm	Code
L1	XXS	55	47	MM-IT-031
L2	XS	60	49	MM-IT-032
L3	S	72	54	MM-IT-033
L4	M	74	57	MM-IT-034
L5	L	76	59	MM-IT-035
L6	XL	78	61	MM-IT-036



Upper




No.	Size	A mm	B mm	Code
U1	XXS	61	47	MM-IT-037
U2	XS	64	51	MM-IT-038
U3	S	66	53	MM-IT-039
U4	M	69	57	MM-IT-040
U5	L	71	60	MM-IT-041
U6	XL	73	62	MM-IT-042

Set	Contents	Code
6 Pcs	L3,L4,L5,U3,U4,U5	MM-IT-043
8 Pcs	L3,L4,L5,L6,U3,U4,U5,U6	MM-IT-044
12 Pcs	L1,L2,L3,L4,L5,L6 U1,U2,U3,U4,U5,U6	MM-IT-045

Solid with Rim




Lower



No.	Size	A mm	B mm	Code
L1	XXS	55	47	MM-IT-046
L2	XS	60	49	MM-IT-047
L3	S	72	54	MM-IT-048
L4	M	74	57	MM-IT-049
L5	L	76	59	MM-IT-050
L6	XL	78	61	MM-IT-051



Upper



No.	Size	A mm	B mm	Code
U1	XXS	61	47	MM-IT-052
U2	XS	64	51	MM-IT-053
U3	S	66	53	MM-IT-054
U4	M	69	57	MM-IT-055
U5	L	71	60	MM-IT-056
U6	XL	73	62	MM-IT-057


Set	Contents	Code
6 Pcs	L3,L4,L5,U3,U4,U5	MM-IT-058
8 Pcs	L3,L4,L5,L6,U3,U4,U5,U6	MM-IT-059
12 Pcs	L1,L2,L3,L4,L5,L6 U1,U2,U3,U4,U5,U6	MM-IT-060

Impression Trays Regular

Perforated with Rim




Lower



No.	Size	A mm	B mm	Code
L1	XL	80	63	MM-IT-061
L2	L	78	59	MM-IT-062
L3	M	76	57	MM-IT-063
L4	NM	72	57	MM-IT-064
L5	S	69	50	MM-IT-065



Upper




No.	Size	A mm	B mm	Code
U1	XL	80	63	MM-IT-066
U2	L	78	59	MM-IT-067
U3	M	76	57	MM-IT-068
U4	NM	72	57	MM-IT-069
U5	S	69	50	MM-IT-070

Set	Contents	Code
10 Pcs	L1,L2,L3,L4,L5, U1,U2,U3,U4,U5	MM-IT-061-070

Solid with Rim




Lower



No.	Size	A mm	B mm	Code
L1	XL	80	63	MM-IT-071
L2	L	78	59	MM-IT-072
L3	M	76	57	MM-IT-073
L4	NM	72	57	MM-IT-074
L5	S	69	50	MM-IT-075



Upper



No.	Size	A mm	B mm	Code
U1	XL	80	63	MM-IT-076
U2	L	78	59	MM-IT-077
U3	M	76	57	MM-IT-078
U4	NM	72	57	MM-IT-079
U5	S	69	50	MM-IT-080


Set	Contents	Code
10 Pcs	L1,L2,L3,L4,L5, U1,U2,U3,U4,U5	MM-IT-071-080

Impression Trays For Children

Perforated with Rim




Lower



No.	Size	A mm	B mm	Code
L1	S	54	47	MM-IT-081
L2	M	62	51	MM-IT-082
L3	L	70	55	MM-IT-083



Upper




No.	Size	A mm	B mm	Code
U1	S	55	43	MM-IT-084
U2	M	60	45	MM-IT-085
U3	L	63	49	MM-IT-086

Set	Contents	Code
6 Pcs	L1,L2,L3,U1,U2,U3	MM-IT-081-086




Lower



No.	Size	A mm	B mm	Code
L1	S	54	47	MM-IT-087
L2	M	62	51	MM-IT-088
L3	L	70	55	MM-IT-089




Upper




No.	Size	A mm	B mm	Code
U1	S	55	43	MM-IT-090
U2	M	60	45	MM-IT-091
U3	L	63	49	MM-IT-092

Set	Contents	Code
6 Pcs	L1,L2,L3,U1,U2,U3	MM-IT-087-092

Impression Trays Regular Ehrlicke Pattern Perforated




Lower



No.	Size	A mm	B mm	Code
UB0	XS	65	46	MM-IT-093
UB1	S	67	50	MM-IT-094
UB2	M	76	54	MM-IT-095
UB3	L	79	56	MM-IT-096
UB4	XL	79	61	MM-IT-097



Upper





No.	Size	A mm	B mm	Code
UB0	XS	59	52	MM-IT-098
UB1	S	66	52	MM-IT-099
UB2	M	71	62	MM-IT-100
UB3	L	71	64	MM-IT-101
UB4	XL	80	70	MM-IT-102

Set	Contents	Code
10 Pcs	LB0,LB1,LB,2LB3,LB4, UB0,UB1,UB2,UB3,UB4	MM-IT-093-102



Impression Trays

Impression Trays Edentulous

Perforated with Rim





Solid with Rim





Impression Trays Anterior Depressed

Perforated with Rim



Solid with Rim



Perforated with Rim



Fig. 30
MM-IT-131
Upper Left / Lower Right



Fig. 31
MM-IT-132
Lower Left / Upper Right



Fig. 32
MM-IT-133
Lower Jaw



Fig. 33
MM-IT-134
Upper Jaw



Fig. 99
MM-IT-135
Lateral

Solid with Rim



Fig. 30
MM-IT-136
Upper Left / Lower Right



Fig. 31
MM-IT-137
Lower Left / Upper Right



Fig. 32
MM-IT-138
Lower Jaw



Fig. 33
MM-IT-139
Upper Jaw



Fig. 99
MM-IT-140
Lateral

Perforated without Rim



Lower Left / Upper Right		
No.	Size	Code
Fig. L2D	S	MM-IT-141
Fig. L4D	M	MM-IT-142
Fig. L6D	L	MM-IT-143



Upper Left / Lower Right		
No.	Size	Code
Fig. L2S	S	MM-IT-147
Fig. L4S	M	MM-IT-148
Fig. L6S	L	MM-IT-149

Perforated with Rim



Lower Left / Upper Right		
No.	Size	Code
Fig. L2D	S	MM-IT-144
Fig. L4D	M	MM-IT-145
Fig. L6D	L	MM-IT-146



Upper Left / Lower Right		
No.	Size	Code
Fig. L2S	S	MM-IT-150
Fig. L4S	M	MM-IT-151
Fig. L6S	L	MM-IT-152

Carvers for silicone impression materials




MM-IT-153




MM-IT-154


Gingival cord packers




MM-IT-155 Serrated



MM-IT-157 Serrated



MM-IT-156 Smooth



MM-IT-158 Serrated

Page

155–164

- SPIRIT LAMPS & JARS
- BOWLS & KIDNEY TRAYS
- INSTRUMENTS TRAYS
- INSTRUMENTS BOXES
- STERILIZING BOXES & INSTRUMENTS TRAY
- MEASURING JUGS, BUCKETS & IRRIGATORS
- NEEDLE CASES, PETRI DISHES & FUNNELS
- INSTRUMENTS TRAYS & ACCESSORIES FOR BURS
- ENDODONTIC TRAYS & INSTRUMENTS

Spirit lamps & Jars




MM-HW-001 Spirit Lamp 120 cc




MM-HW-002 Spirit Lamp 120 cc



MM-HW-003 Spirit Lamp 120 cc



MM-HW-004 Spare wicks Flat 12mm (pack of 5)



MM-HW-005 Spare wicks Round Ø4mm (pack of 5)



Forceps Jars

	dia.	height
MM-HW-006	30	90 mm
MM-HW-007	50	100 mm
MM-HW-008	50	130 mm
MM-HW-009	50	175 mm
MM-HW-010	75	100 mm
MM-HW-011	75	130 mm
MM-HW-012	75	175 mm
MM-HW-013	100	175 mm
MM-HW-014	100	200 mm



Dressing Jars With Lid & Knob

	dia.	height
MM-HW-015	75 mm	60 mm
MM-HW-016	75 mm	75 mm
MM-HW-017	100 mm	60 mm
MM-HW-018	100 mm	100 mm
MM-HW-019	125 mm	80 mm
MM-HW-020	125 mm	125 mm
MM-HW-021	150 mm	150 mm



Dressing Jars / Cotton Containers

MM-HW-022	80 x 60 mm
MM-HW-023	80 x 80 mm
MM-HW-024	100 x 60 mm
MM-HW-025	100 x 80 mm
MM-HW-026	100 x 100 mm
MM-HW-027	120 x 80 mm
MM-HW-028	120 x 100 mm
MM-HW-029	120 x 120 mm
MM-HW-030	150 x 100 mm
MM-HW-031	150 x 150 mm
MM-HW-032	180 x 150 mm
MM-HW-033	180 x 180 mm
MM-HW-034	200 x 200 mm
MM-HW-035	240 x 240 mm



Gallipot	
	vol.
MM-HW-036	50 ml
MM-HW-037	100 ml
MM-HW-038	150 ml
MM-HW-039	225 ml



Stainless Steel Round Bowls			
	dia.	height	vol.
MM-HW-040	Ø 40 mm x 19 mm		0.05L
MM-HW-041	Ø 60 mm x 30 mm		0.07L
MM-HW-042	Ø 80 mm x 40 mm		0.14L
MM-HW-043	Ø 100 mm x 50 mm		0.35L
MM-HW-044	Ø 120 mm x 60 mm		0.45L
MM-HW-045	Ø 150 mm x 70 mm		0.70L
MM-HW-046	Ø 180 mm x 90 mm		1.10L
MM-HW-047	Ø 200 mm x 100 mm		1.60L



Stainless Steel Bowls			
	dia.	height	vol.
MM-HW-048	Ø 40 mm x 19 mm		0.05L
MM-HW-049	Ø 60 mm x 30 mm		0.07L
MM-HW-050	Ø 80 mm x 40 mm		0.14L
MM-HW-051	Ø 100 mm x 50 mm		0.35L
MM-HW-052	Ø 120 mm x 60 mm		0.45L
MM-HW-053	Ø 150 mm x 70 mm		0.70L
MM-HW-054	Ø 180 mm x 90 mm		1.10L
MM-HW-055	Ø 200 mm x 100 mm		1.60L



Wash Basins	
	dia.
MM-HW-056	300 mm
MM-HW-057	320 mm
MM-HW-058	350 mm
MM-HW-059	400 mm



Medicine Cup		
	dia.	height
MM-HW-060	80 x 35 mm	
MM-HW-061	125 x 60 mm	



Medicine Cup (Graduated)		
	Size (mm)	Cap.(ml)
MM-HW-062	2 Oz.	60 ml



Kidney Trays With LID	
	dia.
MM-HW-063	150 mm
MM-HW-064	200 mm
MM-HW-065	250 mm
MM-HW-066	300 mm



Kidney Trays	
	dia.
MM-HW-067	150 mm
MM-HW-068	200 mm
MM-HW-069	250 mm
MM-HW-070	300 mm



Kidney Trays Shallow	
	dia.
MM-HW-071	150 mm
MM-HW-072	200 mm
MM-HW-073	250 mm
MM-HW-074	300 mm

Instruments Trays



Scalers Trays	
MM-HW-075	190 x 100 x 17 mm
MM-HW-076	200 x 100 x 25 mm



Scalers Trays	
MM-HW-077	185 x 80 x 15 mm



Instruments Trays With LID	
MM-HW-078	150 x 100 x 50 mm
MM-HW-079	200 x 100 x 50 mm
MM-HW-080	220 x 120 x 50 mm
MM-HW-081	240 x 140 x 50 mm
MM-HW-082	200 x 150 x 50 mm
MM-HW-083	250 x 200 x 50 mm
MM-HW-084	300 x 250 x 50 mm
MM-HW-085	300 x 200 x 50 mm
MM-HW-086	350 x 250 x 50 mm
MM-HW-087	400 x 300 x 50 mm
MM-HW-088	450 x 300 x 50 mm
MM-HW-089	450 x 350 x 50 mm



Instruments Trays Without LID	
MM-HW-090	150 x 100 x 50 mm
MM-HW-091	200 x 100 x 50 mm
MM-HW-092	220 x 120 x 50 mm
MM-HW-093	240 x 140 x 50 mm
MM-HW-094	200 x 150 x 50 mm
MM-HW-095	250 x 200 x 50 mm
MM-HW-096	300 x 250 x 50 mm
MM-HW-097	300 x 200 x 50 mm
MM-HW-098	350 x 250 x 50 mm
MM-HW-099	400 x 300 x 50 mm
MM-HW-100	450 x 300 x 50 mm
MM-HW-101	450 x 350 x 50 mm



Instruments Trays	
MM-HW-102	210 x 160 x 10 mm
MM-HW-103	240 x 180 x 10 mm
MM-HW-104	310 x 150 x 10 mm
MM-HW-105	310 x 210 x 10 mm
MM-HW-106	350 x 240 x 10 mm
MM-HW-107	400 x 270 x 10 mm
MM-HW-108	420 x 280 x 10 mm



Mayo Trays Perforated	
MM-HW-109	150 x 100 x 50 mm
MM-HW-110	200 x 100 x 50 mm
MM-HW-111	220 x 120 x 50 mm
MM-HW-112	240 x 140 x 50 mm
MM-HW-113	200 x 150 x 50 mm
MM-HW-114	250 x 200 x 50 mm
MM-HW-115	300 x 200 x 50 mm
MM-HW-116	300 x 250 x 50 mm
MM-HW-117	350 x 250 x 50 mm
MM-HW-118	400 x 300 x 50 mm
MM-HW-119	450 x 300 x 50 mm
MM-HW-120	450 x 350 x 50 mm



DIN Baskets perforated sides	
MM-HW-121	240 x 240 x 50 mm
MM-HW-122	240 x 240 x 100 mm
MM-HW-123	480 x 260 x 50 mm
MM-HW-124	520 x 260 x 80 mm
MM-HW-125	520 x 260 x 100 mm



Wiremesh Baskets perforated	
MM-HW-126	240 x 240 x 50 mm
MM-HW-127	260 x 240 x 100 mm
MM-HW-128	480 x 260 x 50 mm
MM-HW-129	520 x 260 x 80 mm
MM-HW-130	520 x 260 x 100 mm

Instruments Boxes



Instrument Boxes With Side Handle	
MM-HW-131	150 x 100 x 50 mm
MM-HW-132	200 x 100 x 50 mm
MM-HW-133	220 x 120 x 50 mm
MM-HW-134	240 x 140 x 50 mm
MM-HW-135	200 x 150 x 50 mm
MM-HW-136	250 x 200 x 50 mm
MM-HW-137	300 x 200 x 50 mm
MM-HW-138	300 x 250 x 50 mm
MM-HW-139	400 x 300 x 50 mm
MM-HW-140	450 x 300 x 50 mm
MM-HW-141	450 x 350 x 50 mm

Instrument Boxes With LID	
MM-HW-142	150 x 100 x 50 mm
MM-HW-143	200 x 100 x 50 mm
MM-HW-144	220 x 120 x 50 mm
MM-HW-145	240 x 140 x 50 mm
MM-HW-146	200 x 150 x 50 mm
MM-HW-147	250 x 200 x 50 mm
MM-HW-148	300 x 200 x 50 mm
MM-HW-149	300 x 250 x 50 mm
MM-HW-150	400 x 300 x 50 mm
MM-HW-151	400 x 300 x 50 mm
MM-HW-152	450 x 350 x 50 mm



Perforated Instruments Boxes	
MM-HW-153	170 x 070 x 35 mm
MM-HW-154	180 x 080 x 40 mm
MM-HW-155	200 x 100 x 40 mm
MM-HW-156	200 x 100 x 50 mm
MM-HW-157	200 x 100 x 60 mm
MM-HW-158	220 x 120 x 60 mm
MM-HW-159	250 x 120 x 60 mm
MM-HW-160	300 x 150 x 60 mm
MM-HW-161	300 x 200 x 80 mm
MM-HW-162	400 x 160 x 80 mm
MM-HW-163	400 x 200 x 100 mm
MM-HW-164	500 x 200 x 120 mm

Perforated Instruments Boxes	
MM-HW-165	170 x 070 x 35 mm
MM-HW-166	180 x 080 x 40 mm
MM-HW-167	200 x 100 x 40 mm
MM-HW-168	200 x 100 x 50 mm
MM-HW-169	200 x 100 x 60 mm
MM-HW-170	220 x 120 x 60 mm
MM-HW-171	250 x 120 x 60 mm
MM-HW-172	300 x 150 x 60 mm
MM-HW-173	300 x 200 x 80 mm
MM-HW-174	400 x 160 x 80 mm
MM-HW-175	400 x 200 x 100 mm
MM-HW-176	500 x 200 x 120 mm



Catheter trays with cover

MM-HW-177	100 x 200 x 50 mm
MM-HW-178	100 x 300 x 50 mm
MM-HW-179	120 x 220 x 50 mm



Instruments Boxes Perforated

MM-HW-180	180 x 80 x 30 mm
MM-HW-181	180 x 80 x 40 mm
MM-HW-182	180 x 80 x 50 mm
MM-HW-183	200 x 90 x 30 mm
MM-HW-184	200 x 90 x 40 mm
MM-HW-185	200 x 90 x 50 mm
MM-HW-186	250 x 100 x 40 mm
MM-HW-187	250 x 100 x 50 mm
MM-HW-188	250 x 100 x 60 mm



Thermometer Jar


	dia. height
MM-HW-189	Ø 25 x 100 mm



Instruments Boxes With Filters

MM-HW-193	170 x 070 x 35 mm
MM-HW-194	180 x 080 x 40 mm
MM-HW-195	200 x 100 x 40 mm
MM-HW-196	200 x 100 x 50 mm
MM-HW-197	200 x 100 x 60 mm
MM-HW-198	220 x 120 x 60 mm
MM-HW-199	250 x 120 x 60 mm
MM-HW-200	300 x 150 x 60 mm
MM-HW-201	300 x 200 x 80 mm
MM-HW-202	400 x 160 x 80 mm
MM-HW-203	400 x 200 x 100 mm
MM-HW-204	500 x 200 x 120 mm

Sterilizing Forceps



Cheatle Forceps

MM-HW-190 27 cm

Collin Forceps

MM-HW-191 20 cm

MM-HW-192 30 cm



Measuring Jugs , Buckets Irrigators



Irrigators	
	vol.
MM-HW-205	1.0L
MM-HW-206	1.5L
MM-HW-207	2.0L



Measuring Jugs	
	vol.
MM-HW-208	250 ml
MM-HW-209	500 ml
MM-HW-210	1000 ml
MM-HW-211	2000 ml



Measuring Jugs	
	vol.
MM-HW-212	250 ml
MM-HW-213	500 ml
MM-HW-214	1000 ml
MM-HW-215	2000 ml



Buckets With LID		Buckets Without LID	
	vol.		vol.
MM-HW-216	8.0L	MM-HW-219	8.0L
MM-HW-217	10.0L	MM-HW-220	10.0L
MM-HW-218	15.0L	MM-HW-221	0L



Sputum Cup	
	dia. height
MM-HW-222	Ø 75 x 80 mm
MM-HW-223	Ø 100 x 50 mm
MM-HW-224	Ø 100 x 80 mm



Ointment Jar		
	dia. height	Vol.
MM-HW-225	Ø 75 x 65 mm	250ml



Needle Case	
MM-HW-226	210 x 125 x 17 mm



Needle Case Sauerbruch	
MM-HW-227	210 x 125 x 17 mm



Needle Case	
MM-HW-228	150 x 90 x 10 mm



Needle Cases Perforated	
MM-HW-229	50 x 25 x 4.5 mm
MM-HW-230	65 x 42 x 8 mm
MM-HW-231	75 x 25 x 4.5 mm



Needle Cases Plain	
MM-HW-232	50 x 25 x 4.5 mm
MM-HW-233	65 x 42 x 8 mm
MM-HW-234	75 x 25 x 4.5 mm



Needle Case Complete	
MM-HW-235	




Funnels	
	dia.
MM-HW-236	100 mm
MM-HW-237	120 mm
MM-HW-238	150 mm
MM-HW-239	200 mm
MM-HW-240	240 mm



Petri Dishes	
	dia. height
MM-HW-241	60 x 20 mm
MM-HW-242	90 x 15 mm
MM-HW-243	100 x 20 mm
MM-HW-244	120 x 20 mm

Instruments Trays & Accessories for Burs



Chico

Stainless Steel for 5 instruments

MM-HW-245 180 x 60 x 22 mm



Medio

Stainless Steel for 10 instruments

MM-HW-246 180 x 125 x 22 mm



Maxi

Stainless Steel for 20 instruments

MM-HW-247 280 x 180 x 25 mm



Maxi Set 1

Stainless Steel rack for 15 instruments

MM-HW-248 288 x 187 x 40 mm



Maxi Set 2

Stainless Steel for 15 instruments


MM-HW-249 288 x 187 x 40 mm



Maxi Set 3

Stainless Steel for 15 instruments

MM-HW-250 288 x 187 x 40 mm



Medio

Stainless Steel Perforated Tray with cover

MM-HW-251 187 x 145 x 25 mm



Medio

Stainless Steel Solid Tray with cover

MM-HW-252 187 x 145 x 25 mm




Burs Disinfecting Box Stainless Steel

MM-HW-253 90 mm x 40 mm



Mesh Burs Holder

MM-HW-254 35 mm x 150 mm



MM-HW-255 Bur Holding Tweezer 16.5 cm



Perforated Instruments Boxes

MM-HW-256	170 x 070 x 35 mm
MM-HW-257	180 x 080 x 40 mm
MM-HW-258	200 x 100 x 40 mm
MM-HW-259	200 x 100 x 50 mm
MM-HW-260	200 x 100 x 60 mm
MM-HW-261	220 x 120 x 60 mm
MM-HW-262	250 x 120 x 60 mm
MM-HW-263	300 x 150 x 60 mm
MM-HW-264	300 x 200 x 80 mm
MM-HW-265	400 x 160 x 80 mm
MM-HW-266	400 x 200 x 100 mm
MM-HW-267	500 x 200 x 120 mm



Stainless Steel endodontic tray with cover

MM-HW-268 250 x 115 x 15 mm



Instrument Box perforated

MM-HW-269 180 x 90 x 30 mm



MM-HW-270 Endometer Laser etched




MM-HW-271 Lester Pumice Cup



MM-HW-272 Endo Box Round



MM-HW-273 Endo Box



MM-HW-274 Endo Box Trio



MM-HW-275 Cotton Holder with Lid



MM-HW-276 Cotton Holder without Lid



MM-HW-277 Cotton Holder Estrella without Lid




Page 165–170

- ORAL SURGICAL BURS
- TRIMMER BURS
- LINDEMANN SIDE CUTTING BURS
- ENT BURS

Toller Round Head Bur

Tungsten Carbide
with Blood Ring




	50mm	70mm
Length	50mm	70mm
Head size	2.3mm	2.3mm
Code No.	MM-BU-001	MM-BU-002

1/1

Sterile Packed in boxes of 10

Toller Taper Fissure Bur

Tungsten Carbide cross cut
with Blood Ring



	50mm	70mm
Length	50mm	70mm
Head size	1.6mm	1.6mm
Code No.	MM-BU-003	MM-BU-004

1/1

Sterile Packed in boxes of 10

Apicectomy Bur

44mm length




Head size	0.8mm
Code No.	MM-BU-005

1/1

Sterile Packed in boxes of 10

Micro Head Apicectomy Bur

10mm length
Latch Fitting



Head size	1.0mm
Code No.	MM-BU-006

1/1







Sterile Packed in boxes of 5








Round Head Bur

Plain cut
44mm length



1/1

						
Head size	1.0mm	1.2mm	1.4mm	1.6mm	1.8mm	2.1mm
Code No.	MM-BU-007	MM-BU-008	MM-BU-009	MM-BU-010	MM-BU-011	MM-BU-012

						
2.3mm	2.5mm	2.7mm	3.1mm	4.0mm	5.0mm	6.0mm
MM-BU-013	MM-BU-014	MM-BU-015	MM-BU-016	MM-BU-017	MM-BU-018	MM-BU-019

Sterile Packed in boxes of 10


Oral Surgery Burs

Flat Fissure Bur

Cross cut
44mm length



1/1



Head size	1mm	1.4mm	1.6mm	1.8mm	2.1mm	2.3mm
Code No.	MM-BU-020	MM-BU-021	MM-BU-022	MM-BU-023	MM-BU-024	MM-BU-025


Sterile Packed in boxes of 10

Tapered Fissure Bur

Cross cut
44mm length



1/1




Head size	1mm	1.4mm	1.6mm	1.8mm	2.1mm	2.3mm
Code No.	MM-BU-026	MM-BU-027	MM-BU-028	MM-BU-029	MM-BU-030	MM-BU-031

Sterile Packed in boxes of 10

Tapered Fissure Bur

Cross cut
50mm length



Head size	1.6mm
Code No.	MM-BU-032

Sterile Packed in boxes of 10


Bud Bur

Cross cut

44mm length 1/1



Head size



5mm

6mm

7mm

44mm length Code No.

MM-BU-033

MM-BU-035

MM-BU-037

70mm length Code No.

MM-BU-034

MM-BU-036

MM-BU-038

44mm Sterile Packed in boxes of 10

70mm Sterile Packed in boxes of 5


Cylinder Bur

Cross cut

44mm length 1/1



Head size



5mm

6mm

7mm

44mm length Code No.

MM-BU-039

MM-BU-041

MM-BU-043

70mm length Code No.

MM-BU-040

MM-BU-042

MM-BU-044

44mm Sterile Packed in boxes of 10

70mm Sterile Packed in boxes of 5


Pear Bur

Cross cut

44mm length 1/1



Head size



5mm

6mm

7mm

44mm length Code No.

MM-BU-045

MM-BU-047

MM-BU-049

70mm length Code No.

MM-BU-046

MM-BU-048








MM-BU-050

44mm Sterile Packed in boxes of 10

70mm Sterile Packed in boxes of 5

Lindemann Side Cutting Burs








Lindemann Side Cutting Burs

						
	Length 44mm		65mm		75mm	
	Code No. MM-BU-051		MM-BU-052		MM-BU-053	
		1/1		1/1		1/1

Sterile Packed in boxes of 5

Bone Trephines

Bone Trephines

						
		Length 44mm	44mm	44mm		44mm
		Head size 1.8mm	2.3mm	2.7mm		5mm
		Code No. MM-BU-054	MM-BU-055	MM-BU-056		MM-BU-057
1/1					1/1	

Sterile Packed in boxes of 5








Sterile Packed in boxes of 5









ENT Cutting Bur

70mm length



1/1

							
Head size	0.6mm	0.7mm	0.8mm	1mm	1.4mm	1.8mm	2.3mm
Code No.	MM-BU-058	MM-BU-059	MM-BU-060	MM-BU-061	MM-BU-062	MM-BU-063	MM-BU-064

							
2.7mm	3.1mm	3.5mm	4mm	4.5mm	5mm	6mm	7mm
MM-BU-065	MM-BU-066	MM-BU-067	MM-BU-068	MM-BU-069	MM-BU-070	MM-BU-071	MM-BU-072

Sterile Packed in boxes of 5

ENT Allport Cutting Bur








70mm length









6 Blades



1/1

Bur

							
Head size	0.6mm	0.7mm	0.8mm	1mm	1.4mm	1.8mm	2.3mm
Code No.	MM-BU-073	MM-BU-074	MM-BU-075	MM-BU-076	MM-BU-077	MM-BU-078	MM-BU-079

							
2.7mm	3.1mm	3.5mm	4mm	4.5mm	5mm	6mm	7mm
MM-BU-080	MM-BU-081	MM-BU-082	MM-BU-083	MM-BU-084	MM-BU-085	MM-BU-086	MM-BU-087

Sterile Packed in boxes of 5








ENT Burs









ENT Polishing Bur

70mm length



1/1

							
Head size	0.6mm	0.7mm	0.8mm	1mm	1.4mm	1.8mm	2.3mm
Code No.	MM-BU-088	MM-BU-089	MM-BU-090	MM-BU-091	MM-BU-092	MM-BU-093	MM-BU-094

							
2.7mm	3.1mm	3.5mm	4mm	4.5mm	5mm	6mm	7mm
MM-BU-095	MM-BU-096	MM-BU-097	MM-BU-098	MM-BU-099	MM-BU-100	MM-BU-101	MM-BU-102








Sterile Packed in boxes of 5









ENT Diamond Bur

70mm length



1/1

							
Head size	0.6mm	0.7mm	0.8mm	1mm	1.4mm	1.8mm	2.3mm
Code No.	MM-BU-103	MM-BU-104	MM-BU-105	MM-BU-106	MM-BU-107	MM-BU-108	MM-BU-109

							
2.7mm	3.1mm	3.5mm	4mm	4.5mm	5mm	6mm	7mm
MM-BU-110	MM-BU-111	MM-BU-112	MM-BU-113	MM-BU-114	MM-BU-115	MM-BU-116	MM-BU-117

Sterile Packed in boxes of 5

Registered Office:

MAYOMED LIMITED

8b, Kelvin House, Kelvin Way
Carwley, West Sussex, RH10 9WE
United Kingdom

T: +44 (0) 1293 344 447

F: +44 (0) 1293 344 448

M: +44 (0) 7951 900 642

E: info@mayo-med.com

Web: www.mayo-med.com

Company Registration No. 11578206
Registered in England & Wales



ISO 9001:2015
ISO 13485:2016